



 **Bledsoe**

BRACING FOR THE FUTURE





## TABLE OF CONTENTS

---



- 6** OA  
THRUSTER  
OA AIR  
ALIGNER  
AXIOM OA  
Z-12 OA
- 10** FUNCTIONAL  
Z-12  
AXIOM-D  
AXIOM  
ULTIMATE DYNAMIC  
ULTIMATE SPORT  
CROSSOVER  
KNEE SUPPORT
- 16** PATELLOFEMORAL  
20.50  
REINFORCED PATELLA STABILIZER  
LATERAL PATELLA STABILIZER  
TENDON STRAP
- 20** POST-OP  
EXTENDER PLUS  
MERIT OR  
ORIGINAL  
LEVER LOCK  
CROSSOVER ROM  
THREE PANEL KNEE IMMOBILIZER  
SINGLE PANEL KNEE IMMOBILIZER  
PHILIPPON HIP  
SIMPLE PELVIC
- 26** FOOT AND ANKLE  
AXIOM ANKLE  
ULTIMATE ANKLE  
SIDEKICK  
WRAPTOR ANKLE  
M-BRACE
- 30** WALKER BOOTS  
EVERYDAY BOOTS  
PEDIATRIC BOOT  
LONG-TERM BOOTS  
DIABETIC BOOTS  
SPECIALTY BOOTS  
PFS  
SPLINTS  
SILICONE HEEL CUPS
- 42** UPPER EXTREMITY  
ARC / ARC-XR  
ARC WITH PILLOW  
CLINIC ARM SLING  
CLINIC SHOULDER IMMOBILIZER  
UNIVERSAL ARM SLING  
SLING AND SWATHE IMMOBILIZER  
CLAVICLE SUPPORT
- 46** ELBOW, WRIST AND HAND  
TELESCOPING ELBOW  
TENDON STRAP  
WRIST LACER  
WRIST LACER WITH THUMB SPICA  
WRIST BRACE  
WRIST BRACE WITH THUMB SPICA  
TRIGGER AND KEEPER THUMB SUPPORT
- 50** SPINE  
VOLARE WITH CHAIRBACK  
VOLARE LS—XT  
VOLARE LS  
VOLARE LITE  
VOLARE MODULAR SYSTEM
- 54** COLD THERAPY  
COLD CONTROL  
UNIVERSAL PAD  
SHOULDER PAD  
HAND, WRIST AND FINGERS PAD  
ANKLE PAD  
5X10 RECTANGLE PAD  
10X14 RECTANGLE PAD  
SHOULDER WRAP
- 58** INDEX
- 62** WARRANTY/ORDERING INFO






## FROM GREG NELSON

---

**D**ear valued customer:

In this day of economic uncertainty and our country's time of critical change, sometimes we forget to count our blessings and have some fun! I've enjoyed adding a little "color" to our products in this catalog; I hope you enjoy this year's cover and theme: *Bracing for the Future!*

At Bledsoe, we work hard to engineer and manufacture the best products. We realize that somebody very special will be using our product to treat their injury. We also take our responsibility to the medical community extremely seriously; knowing that you have chosen our product to treat your patients.



From all of our Bledsoe team, I want to say thanks for choosing Bledsoe and I promise—we will not let you down!

Sincerely and Best Wishes,


---

Chairman, Chief Executive Officer  
Bledsoe Brace Systems




## COLOR SELECTION SYSTEM


Select one or two colors for your brace. Fill in the code for your color selection on the Custom Brace Order Form provided by your fitter. For a single color, use the same code twice. Color options are only available on custom knee braces.




Signal Yellow  
1003




Dahlia Yellow  
1037




Hot Pink  
4010




Traffic Blue  
5017




Tan  
1001




Pure Orange  
2004




Gentian Blue  
5002




Water Blue  
5021




Canary  
1016




Brick  
3011




Aquablue  
5009




Gray Blue  
5024




Zinc Yellow  
1018




Pure Violet  
4008



Sky Blue  
5015




Violet Blue  
5000




Racing Red  
3020









Moss Green  
6005




Jet Black  
9005



Sparkle Green  
50024




Lawn  
6010




Pure White  
9016




Stardust  
Blue/Purple  
40001




Yellow Green  
6018




Dark Gray  
7016




Sparkle  
Sky Blue  
40074




Light Gray  
7032




Dusty Gray  
7043



Sparkle Silver  
90450



Sparkle Black  
88420



Sparkle Red  
30028





DA




THRUSTER RLF

DA AIR

ALIGNER ESE

AXIOM DA

Z-12 DA





## BLED SOE PRINCIPLE:

A BETTER OA SOLUTION

Our proven technology works.

It relieves pain.

"The Thruster remains the best performing OA brace on the market today."

— *In Vivo Three Dimensional Determination of OA Brace Effectiveness: A Multiple Brace Analysis,* Richard D. Komistek, Ph.D., Director, Center for Musculoskeletal Research, *Journal of Bone and Joint Surgery* 2005; 87; 114-119.



VARUS ALIGNMENT  
OF THE KNEE



FORCES APPLIED BY  
THRUSTER TECHNOLOGY

## THRUSTER RLF

- > Dynamic control provides the greatest relief in full extension and relaxes in flexion
- > Eliminates pressure points and increases patient compliance
- > The only single upright brace proven to provide pain relief\*
- > Patient adjustable dial

### Indications:

For medial or lateral Unicompartmental Osteoarthritis.

- Custom: Fabricated from FitKit 360° measurements, laser scan, or cast mold.
- Off-the-Shelf: Measure 6" above and 6" below mid-patella.

Thigh Circumference	Calf Circumference				Medial		Lateral	
	STD and Short	Athletic			Left	Right	Left	Right
15.0" - 17.5" (38.1 cm - 44.5 cm)	12.0" - 13.5" (30.5 cm - 34.3 cm)	10.5" - 12.0" (26.7 cm - 30.5 cm)	XS	STD	TA009101--B	TA009201--B	TA008101--B	TA008201--B
				Short	TA019101--B	TA019201--B	TA018101--B	TA018201--B
				Athletic	TA029101--B	TA029201--B	TA028101--B	TA028201--B
17.5" - 20.0" (44.5 cm - 50.8 cm)	13.5" - 15.0" (34.3 cm - 38.1 cm)	12.0" - 13.5" (30.5 cm - 34.3 cm)	S	STD	TA009103--B	TA009203--B	TA008103--B	TA008203--B
				Short	TA019103--B	TA019203--B	TA018103--B	TA018203--B
				Athletic	TA029103--B	TA029203--B	TA028103--B	TA028203--B
20.0" - 22.5" (50.8 cm - 57.2 cm)	15.0" - 16.5" (38.1 cm - 41.9 cm)	13.5" - 15.0" (34.3 cm - 38.1 cm)	M	STD	TA009105--B	TA009205--B	TA008105--B	TA008205--B
				Short	TA019105--B	TA019205--B	TA018105--B	TA018205--B
				Athletic	TA029105--B	TA029205--B	TA028105--B	TA028205--B
22.5" - 25.0" (57.2 cm - 63.5 cm)	16.5" - 18.0" (41.9 cm - 45.7 cm)	15.0" - 16.5" (38.1 cm - 41.9 cm)	L	STD	TA009107--B	TA009207--B	TA008107--B	TA008207--B
				Short	TA019107--B	TA019207--B	TA018107--B	TA018207--B
				Athletic	TA029107--B	TA029207--B	TA028107--B	TA028207--B
25.0" - 27.5" (63.5 cm - 69.9 cm)	18.0" - 19.5" (45.7 cm - 49.5 cm)	16.5" - 18.0" (41.9 cm - 45.7 cm)	XL	STD	TA009109--B	TA009209--B	TA008109--B	TA008209--B
				Short	TA019109--B	TA019209--B	TA018109--B	TA018209--B
				Athletic	TA029109--B	TA029209--B	TA028109--B	TA028209--B
27.5" - 29.0" (69.9 cm - 73.7 cm)	19.5" - 21.0" (49.5 cm - 53.3 cm)	18.0" - 19.5" (45.7 cm - 49.5 cm)	2XL	STD	TA009111--B	TA009211--B	TA008111--B	TA008211--B
				Short	TA019111--B	TA019211--B	TA018111--B	TA018211--B
				Athletic	TA029111--B	TA029211--B	TA028111--B	TA028211--B
29.0" - 31.0" (73.7 cm - 78.7 cm)	21.0" - 22.5" (53.3 cm - 57.2 cm)	19.5" - 21.0" (49.5 cm - 53.3 cm)	3XL	STD	TA009113--B	TA009213--B	TA008113--B	TA008213--B
				Short	TA019113--B	TA019213--B	TA018113--B	TA018213--B
				Athletic	TA029113--B	TA029213--B	TA028113--B	TA028213--B

Suggested HCFA Code: L1844 (Custom), L1843 (OTS)

\*Komistek, Richard D. Ph.D., Director, Center for Musculoskeletal Research, Nadaud, Matthew C. MD; Dennis, Douglas, A. MD; Anderle, Matthew R.; Kubo, Matthew B. "In Vivo Three Dimensional Determination of OA Brace Effectiveness: A Multiple Brace Analysis." *American Academy of Orthopedic Surgeons*, February 2005.



Bledsoe  
**Thruster**  
RLF

2009™



## OA AIR

- > Applies the principles of the Thruster in a lightweight, low-profile, shorter brace
- > New unique pneumatic adjustable pad gently increases correction for desired pain relief and comfort
- > Magnesium Material
- > Ideal for patients with mild to moderate OA

### Indications:

Medial or lateral Unicompartmental Osteoarthritis.

- Custom: Fabricated from FitKit 360° measurements, laser scan, or cast mold.
- Off-the-Shelf: Measure 6" above and 6" below mid-patella.

OTS Available  
Fall 2009

Thigh Circumference	Calf Circumference		Medial		Lateral	
			Left	Right	Left	Right
15.0" - 17.5" (38.1 cm - 44.5 cm)	12.0" - 13.5" (30.5 cm - 34.3 cm)	XS	LT009101	LT009201	LT008101	LT008201
17.5" - 20.0" (44.5 cm - 50.8 cm)	13.5" - 15.0" (34.3 cm - 38.1 cm)	S	LT009103	LT009203	LT008103	LT008203
20.0" - 22.5" (50.8 cm - 57.2 cm)	15.0" - 16.5" (38.1 cm - 41.9 cm)	M	LT009105	LT009205	LT008105	LT008205
22.5" - 25.0" (57.2 cm - 63.5 cm)	16.5" - 18.0" (41.9 cm - 45.7 cm)	L	LT009107	LT009207	LT008107	LT008207
25.0" - 27.5" (63.5 cm - 69.9 cm)	18.0" - 19.5" (45.7 cm - 49.5 cm)	XL	LT009109	LT009209	LT008109	LT008209
27.5" - 29.0" (63.5 cm - 69.9 cm)	19.5" - 21.0" (45.7 cm - 49.5 cm)	2XL	LT009111	LT009211	LT008111	LT008211

Suggested HCFA Code: L1844 (Custom), L1843 (OTS)

OA AIR



## ALIGNER ESE

- > Applies the principles of the Thruster with a lower profile hinge design
- > Ideal for patients with mild to moderate OA

### Indications:

Medial or lateral Unicompartmental Osteoarthritis.

- Off-the-Shelf: Measure 6" above and 6" below mid-patella.

Thigh Circumference	Calf Circumference		Medial		Lateral	
			Left	Right	Left	Right
15.0" - 17.5" (38.1 cm - 44.5 cm)	12.0" - 13.5" (30.5 cm - 34.3 cm)	XS	AR079101	AR079201	AR078101	AR078201
17.5" - 20.0" (44.5 cm - 50.8 cm)	13.5" - 15.0" (34.3 cm - 38.1 cm)	S	AR079103	AR079203	AR078103	AR078203
20.0" - 22.5" (50.8 cm - 57.2 cm)	15.0" - 16.5" (38.1 cm - 41.9 cm)	M	AR079105	AR079205	AR078105	AR078205
22.5" - 25.0" (57.2 cm - 63.5 cm)	16.5" - 18.0" (41.9 cm - 45.7 cm)	L	AR079107	AR079207	AR078107	AR078207
25.0" - 27.5" (63.5 cm - 69.9 cm)	18.0" - 19.5" (45.7 cm - 49.5 cm)	XL	AR079109	AR079209	AR078109	AR078209

Suggested HCFA Code: L1843

Bledsoe  
**ALIGNER ESE**





## AXIOM OA

- > Provides mild to moderate OA pain relief while stabilizing the knee and protecting the ACL, PCL, MCL and LCL after injury or repair
- > Up to 8° of varus or valgus correction
- > Optional dynamic version adds posterior force to control tibial translation for patients with ACL or PCL deficient knees
- > Magnesium or Aluminum Material

### Indications:

Mild to moderate medial or lateral Unicompartmental Osteoarthritis while providing support of the ACL, PCL, MCL or LCL and menisci.

- Custom: Fabricated from FitKit 360° measurements, laser scan, or cast mold.

Suggested HCFA Code: L1846

AXIOM OA



## Z-12 OA

- > Provides mild to moderate OA pain relief while stabilizing the knee and protecting the ACL, PCL, MCL and LCL after injury or repair
- > Standard 13" length and EXT 15" length available
- > Magnesium Material

### Indications:

Mild to moderate medial or lateral Unicompartmental Osteoarthritis while providing support of the ACL, PCL, MCL or LCL and menisci.

- Custom: Fabricated from FitKit 360° measurements, laser scan, or cast mold.
- Off-the-Shelf: Measure 6" above and 6" below mid-patella.

Z<sup>12</sup>



Thigh Circumference		Calf Circumference		Medial				Lateral	
	STD and Short	Athletic		Left	Right	Left	Right		
15.0" - 17.5" (38.1 cm - 44.5 cm)	12.0" - 13.5" (30.5 cm - 34.3 cm)	10.5" - 12.0" (26.7 cm - 30.5 cm)	XS	STD 13" EXT 15"	AZ294101 AZ293101	AZ294201 AZ293201	AZ284101 AZ283101	AZ284201 AZ283201	
			Athletic	STD 13" EXT 15"	AZ299101 AZ298101	AZ299201 AZ298201	AZ289101 AZ288101	AZ289201 AZ288201	
17.5" - 20.0" (44.5 cm - 50.8 cm)	13.5" - 15.0" (34.3 cm - 38.1 cm)	12.0" - 13.5" (30.5 cm - 34.3 cm)	S	STD 13" EXT 15"	AZ294103 AZ293103	AZ294203 AZ293203	AZ284103 AZ283103	AZ284203 AZ283203	
			Athletic	STD 13" EXT 15"	AZ299103 AZ298103	AZ299203 AZ298203	AZ289103 AZ288103	AZ289203 AZ288203	
20.0" - 22.5" (50.8 cm - 57.2 cm)	15.0" - 16.5" (38.1 cm - 41.9 cm)	13.5" - 15.0" (34.3 cm - 38.1 cm)	M	STD 13" EXT 15"	AZ294105 AZ293105	AZ294205 AZ293205	AZ284105 AZ283105	AZ284205 AZ283205	
			Athletic	STD 13" EXT 15"	AZ299105 AZ298105	AZ299205 AZ298205	AZ289105 AZ288105	AZ289205 AZ288205	
22.5" - 25.0" (57.2 cm - 63.5 cm)	16.5" - 18.0" (41.9 cm - 45.7 cm)	15.0" - 16.5" (38.1 cm - 41.9 cm)	L	STD 13" EXT 15"	AZ294107 AZ293107	AZ294207 AZ293207	AZ284107 AZ283107	AZ284207 AZ283207	
			Athletic	STD 13" EXT 15"	AZ299107 AZ298107	AZ299207 AZ298207	AZ289107 AZ288107	AZ289207 AZ288207	
25.0" - 27.5" (63.5 cm - 69.9 cm)	18.0" - 19.5" (45.7 cm - 49.5 cm)	16.5" - 18.0" (41.9 cm - 45.7 cm)	XL	STD 13" EXT 15"	AZ294109 AZ293109	AZ294209 AZ293209	AZ284109 AZ283109	AZ284209 AZ283209	
			Athletic	STD 13" EXT 15"	AZ299109 AZ298109	AZ299209 AZ298209	AZ289109 AZ288109	AZ289209 AZ288209	
27.5" - 29.0" (69.9 cm - 73.7 cm)	19.5" - 21.0" (49.5 cm - 53.3 cm)	18.0" - 19.5" (45.7 cm - 49.5 cm)	2XL	STD 13" EXT 15"	AZ294111 AZ293111	AZ294211 AZ293211	AZ284111 AZ283111	AZ284211 AZ283211	
			Athletic	STD 13" EXT 15"	AZ299111 AZ298111	AZ299211 AZ298211	AZ289111 AZ288111	AZ289211 AZ288211	

For PCL/Combined Instability version, indicate condition when ordering.

Suggested HCFA Code: L1846 (Custom), L1845 (OTS)



# FUNCTIONAL



Z-12

AXIOM-D

AXIOM

ULTIMATE DYNAMIC

ULTIMATE SPORT

CROSSOVER

KNEE SUPPORT





## BLED SOE PRINCIPLE:

USE THE BEST MATERIALS

When engineers in the aerospace and automobile industries need an alloy to reduce weight yet maintain strength under duress, they turn to Z-12 Magnesium...

SO DID WE!

### STRONG.



Aerospace

### POWERFUL.



Automotive

### EXTREMELY LIGHTWEIGHT.



Bicycles

## Z-12

- > Most streamlined ligament brace available—13 inch length
- > Strong and durable Z-12 Magnesium
- > Quick release buckles for easy application
- > EXT 15 inch length also available

### Indications:

Provides stability and support following injury or reconstruction of the ACL, PCL, MCL or LCL and menisci.

Z<sup>12</sup>

- Custom: Fabricated from FitKit 360° measurements, laser scan, or cast mold.
- Off-the-Shelf: Measure 6" above and 6" below mid-patella.

Thigh Circumference	Calf Circumference				Standard		Athletic	
	STD	Athletic			Left	Right	Left	Right
13.5" - 16.0" (34.3 cm - 40.6 cm)	12.5" - 14.0" (31.7 cm - 35.6 cm)	11.0" - 12.5" (27.9 cm - 31.7 cm)	XS	STD	AZ114101	AZ114201	AZ113101	AZ113201
16.0" - 18.75" (40.6 cm - 47.6 cm)	14.0" - 15.5" (35.6 cm - 39.4 cm)	12.5" - 14.0" (31.7 cm - 35.6 cm)	S	STD	AZ114103	AZ114203	AZ113103	AZ113203
18.75" - 21.5" (47.6 cm - 54.6 cm)	15.5" - 17.0" (39.4 cm - 43.2 cm)	14.0" - 15.5" (35.6 cm - 39.4 cm)	M	STD	AZ114105	AZ114205	AZ113105	AZ113205
21.5" - 24.25" (54.6 cm - 61.6 cm)	17.0" - 18.5" (43.2 cm - 47.0 cm)	15.5" - 17.0" (39.4 cm - 43.2 cm)	L	STD	AZ114107	AZ114207	AZ113107	AZ113207
24.25" - 27.0" (61.6 cm - 68.6 cm)	18.5" - 20.0" (47.0 cm - 50.8 cm)	17.0" - 18.5" (43.2 cm - 47.0 cm)	XL	STD	AZ114109	AZ114209	AZ113109	AZ113209
27.0" - 29.5" (68.6 cm - 74.9 cm)	20.0" - 21.375" (50.8 cm - 54.3 cm)	18.5" - 20.0" (47.0 cm - 50.8 cm)	2XL	STD	AZ114111	AZ114211	AZ113111	AZ113211
				EXT	AZ119111	AZ119211	AZ118111	AZ118211

For PCL/Combined Instability version, indicate condition when ordering.

Suggested HCFA Code: L1846 (Custom), L1845 (OTS)





## BLED SOE PRINCIPLE:


### DYNAMIC CONTROL

**A**nterior drawer leaves patients with a dreadful feeling of instability. By controlling tibial translation, the Axiom-D helps patients regain confidence to golf, ski and enjoy other lifetime sports.

Our proprietary Dynamic Tibial Mechanism moves posterior as the leg extends controlling tibial subluxation and that feeling of instability.



NO FORCE IN FLEXION



HIGH FORCE IN EXTENSION

## AXIOM-D

- > Dynamic control of the tibia
- > Unique strapping system maintains brace position
- > Maximum strength with our spring steel reinforced Aluminum version
- > Lightweight Magnesium material available Summer 2009

### Indications:

Support for ACL or PCL deficient patients suffering instability, symptoms of giving way, poor quadriceps or hamstring strength or seeking an early return to activities after injury or repair.

- Custom: Fabricated from FitKit 360° measurements, laser scan, or cast mold.
- Off-the-Shelf: Measure 6" above and 6" below mid-patella.

Thigh Circumference	Calf Circumference				Aluminium		Magnesium	
	STD	Athletic			Left	Right	Left	Right
13.5" - 16.0" (34.3 cm - 40.6 cm)	12.5" - 14.0" (31.7 cm - 35.6 cm)	11.0" - 12.5" (27.9 cm - 31.7 cm)	XS	STD Athletic	PK118101 PK113101	PK118201 PK113201	MG118101 MG113101	MG118201 MG113201
16.0" - 18.75" (40.6 cm - 47.6 cm)	14.0" - 15.5" (35.6 cm - 39.4 cm)	12.5" - 14.0" (31.7 cm - 35.6 cm)	S	STD Athletic	PK118103 PK113103	PK118203 PK113203	MG118103 MG113103	MG118203 MG113203
18.75" - 21.5" (47.6 cm - 54.6 cm)	15.5" - 17.0" (39.4 cm - 43.2 cm)	14.0" - 15.5" (35.6 cm - 39.4 cm)	M	STD Athletic	PK118105 PK113105	PK118205 PK113205	MG118105 MG113105	MG118205 MG113205
21.5" - 24.25" (54.6 cm - 61.6 cm)	17.0" - 18.5" (43.2 cm - 47.0 cm)	15.5" - 17.0" (39.4 cm - 43.2 cm)	L	STD Athletic	PK118107 PK113107	PK118207 PK113207	MG118107 MG113107	MG118207 MG113207
24.25" - 27.0" (61.6 cm - 68.6 cm)	18.5" - 20.0" (47.0 cm - 50.8 cm)	17.0" - 18.5" (43.2 cm - 47.0 cm)	XL	STD Athletic	PK118109 PK113109	PK118209 PK113209	MG118109 MG113109	MG118209 MG113209
27.0" - 29.5" (68.6 cm - 74.9 cm)	20.0" - 21.375" (50.8 cm - 54.3 cm)	18.5" - 20.0" (47.0 cm - 50.8 cm)	2XL	STD Athletic	PK118111 PK113111	PK118211 PK113211	MG118111 MG113111	MG118211 MG113211
29.5" - 31.0" (75.0 cm - 79.0 cm)	21.375" - 22.75" (54.3 cm - 57.8 cm)	20.0" - 21.375" (50.8 cm - 54.3 cm)	3XL	STD Athletic	PK118113 PK113113	PK118213 PK113213	MG118113 MG113113	MG118213 MG113213

For PCL/Combined Instability version, indicate condition when ordering.

Suggested HCFA Code: L1846 (Custom), L1845 (OTS)

AXIOM-D



2009™  
**Bledsoe**



# AXIOM

- > Same features and benefits as our AXIOM-D, but for patients who do not need additional dynamic control
- > Available in Magnesium and Aluminum

## AXIOM

### Indications:

Provides stability and support following injury or reconstruction of the ACL, PCL, MCL, LCL, and menisci. Also recommended for prophylactic use.

● Custom: Fabricated from FitKit 360° measurements, laser scan, or cast mold.

● Off-the-Shelf: Measure 6" above and 6" below mid-patella.

Thigh Circumference	Calf Circumference				Aluminium		Magnesium	
	STD	Athletic			Left	Right	Left	Right
13.5" - 16.0" (34.3 cm - 40.6 cm)	12.5" - 14.0" (31.7 cm - 35.6 cm)	11.0" - 12.5" (27.9 cm - 31.7 cm)	XS	STD Athletic	PK114101 PK119101	PK114201 PK119201	MG114101 MG119101	MG114201 MG119201
16.0" - 18.75" (40.6 cm - 47.6 cm)	14.0" - 15.5" (35.6 cm - 39.4 cm)	12.5" - 14.0" (31.7 cm - 35.6 cm)	S	STD Athletic	PK114103 PK119103	PK114203 PK119203	MG114103 MG119103	MG114203 MG119203
18.75" - 21.5" (47.6 cm - 54.6 cm)	15.5" - 17.0" (39.4 cm - 43.2 cm)	14.0" - 15.5" (35.6 cm - 39.4 cm)	M	STD Athletic	PK114105 PK119105	PK114205 PK119205	MG114105 MG119105	MG114205 MG119205
21.5" - 24.25" (54.6 cm - 61.6 cm)	17.0" - 18.5" (43.2 cm - 47.0 cm)	15.5" - 17.0" (39.4 cm - 43.2 cm)	L	STD Athletic	PK114107 PK119107	PK114207 PK119207	MG114107 MG119107	MG114207 MG119207
24.25" - 27.0" (61.6 cm - 68.6 cm)	18.5" - 20.0" (47.0 cm - 50.8 cm)	17.0" - 18.5" (43.2 cm - 47.0 cm)	XL	STD Athletic	PK114109 PK119109	PK114209 PK119209	MG114109 MG119109	MG114209 MG119209
27.0" - 29.5" (68.6 cm - 74.9 cm)	20.0" - 21.375" (50.8 cm - 54.3 cm)	18.5" - 20.0" (47.0 cm - 50.8 cm)	2XL	STD Athletic	PK114111 PK119111	PK114211 PK119211	MG114111 MG119111	MG114211 MG119211
29.5" - 31" (74.9 cm - 78.7 cm)	21.375" - 22.75" (54.3 cm - 57.8 cm)	20.0" - 21.375" (50.8 cm - 54.3 cm)	3XL	STD Athletic	PK114113 PK119113	PK114213 PK119213	----- -----	----- -----

For PCL/Combined Instability version, indicate condition when ordering.

Suggested HCFA Code: L1846 (Custom), L1845 (OTS)





## ULTIMATE DYNAMIC

- > Dynamic control of the tibia
- > Easily formed and adapted in the field
- > Squared frame provides better soft tissue control
- > Maximum strength with our spring steel reinforced Aluminum

### Indications:

Support for ACL or PCL deficient patients suffering instability, symptoms of giving way, poor quadriceps or hamstring strength or seeking an early return to activities after injury or repair.

● Custom: Fabricated from FitKit 360° measurements, laser scan, or cast mold.

● Off-the-Shelf: Measure 6" above and 6" below mid-patella.

Thigh Circumference	Calf Circumference		Everyday Pads		High Activity Pads	
			Left	Right	Left	Right
13.5" - 16.0" (34.3 cm - 40.6 cm)	12.5" - 14.0" (31.7 cm - 35.6 cm)	XS	GK060101	GK060201	JK220101	JK220201
16.0" - 18.75" (40.6 cm - 47.6 cm)	14.0" - 15.5" (35.6 cm - 39.4 cm)	S	GK060103	GK060203	JK220103	JK220203
18.75" - 21.5" (47.6 cm - 54.6 cm)	15.5" - 17.0" (39.4 cm - 43.2 cm)	M	GK060105	GK060205	JK220105	JK220205
21.5" - 24.25" (54.6 cm - 61.6 cm)	17.0" - 18.5" (43.2 cm - 47.0 cm)	L	GK060107	GK060207	JK220107	JK220207
24.25" - 27.0" (61.6 cm - 68.6 cm)	18.5" - 20.0" (47.0 cm - 50.8 cm)	XL	GK060109	GK060209	JK220109	JK220209
27.0" - 29.5" (68.6 cm - 74.9 cm)	20.0" - 21.375" (50.8 cm - 54.3 cm)	2XL	GK060111	GK060211	JK220111	JK220211
29.5" - 31.0" (75.0 cm - 79.0 cm)	21.375" - 22.75" (54.3 cm - 57.8 cm)	3XL	GK060113	GK060213	JK220113	JK220213

Suggested HCFA Code: L1846 (Custom), L1845 (OTS)



## ULTIMATE SPORT

- > Same features and benefits as our Ultimate Dynamic, but for patients who do not need additional dynamic control
- > Maximum strength with our spring steel reinforced Hi-Impact version

### Indications:

Provides stability and support following injury or reconstruction of the ACL, PCL, MCL, LCL, and menisci.

● Custom: Fabricated from FitKit 360° measurements, laser scan, or cast mold.

● Off-the-Shelf: Measure 6" above and 6" below mid-patella.

Thigh Circumference	Calf Circumference			Everyday Pads		High Activity Pads	
				Left	Right	Left	Right
13.5" - 16.0" (34.3 cm - 40.6 cm)	12.5" - 14.0" (31.7 cm - 35.6 cm)	XS	STD Hi-Impact	GK040101 GK070101	GK040201 GK070201	JK020101 JK120101	JK020201 JK120201
16.0" - 18.75" (40.6 cm - 47.6 cm)	14.0" - 15.5" (35.6 cm - 39.4 cm)	S	STD Hi-Impact	GK040103 GK070103	GK040203 GK070203	JK020103 JK120103	JK020203 JK120203
18.75" - 21.5" (47.6 cm - 54.6 cm)	15.5" - 17.0" (39.4 cm - 43.2 cm)	M	STD Hi-Impact	GK040105 GK070105	GK040205 GK070205	JK020105 JK120105	JK020205 JK120205
21.5" - 24.25" (54.6 cm - 61.6 cm)	17.0" - 18.5" (43.2 cm - 47.0 cm)	L	STD Hi-Impact	GK040107 GK070107	GK040207 GK070207	JK020107 JK120107	JK020207 JK120207
24.25" - 27.0" (61.6 cm - 68.6 cm)	18.5" - 20.0" (47.0 cm - 50.8 cm)	XL	STD Hi-Impact	GK040109 GK070109	GK040209 GK070209	JK020109 JK120109	JK020209 JK120209
27.0" - 29.5" (68.6 cm - 74.9 cm)	20.0" - 21.375" (50.8 cm - 54.3 cm)	2XL	STD Hi-Impact	GK040111 GK070111	GK040211 GK070211	JK020111 JK120111	JK020211 JK120211
29.5" - 31.0" (75.0 cm - 79.0 cm)	21.375" - 22.75" (54.3 cm - 57.8 cm)	3XL	STD Hi-Impact	GK040113 GK070113	GK040213 GK070213	JK020113 JK120113	JK020213 JK120213

Suggested HCFA Code: L1846 (Custom), L1845 (OTS)





## CROSSOVER

- > Cool breathable Spacer Mesh provides support, compression and warmth to the knee
- > Easily adjustable hinges for ideal anatomical placement
- > Available in short and standard lengths
- > Also available in Neoprene

### Indications:

For mild ACL, PCL, and MCL sprains, menisci and patellar retinaculum injuries, and mild instabilities.

### Crossover - Spacer Mesh

Thigh Circumference			Pull-Up		Front Closure	
			w/o Buttress	w/Buttress	w/o Buttress	w/Buttress
15.0" - 17.5" (38.1 cm - 44.5 cm)	XS	STD	RK061001	RK062001	RK060001	
		Short	RK064001	RK065001	RK063001	
17.5" - 20.0" (44.5 cm - 50.8 cm)	S	STD	RK061003	RK062003	RK060003	
		Short	RK064003	RK065003	RK063003	
20.0" - 22.5" (50.8 cm - 57.2 cm)	M	STD	RK061005	RK062005	RK060005	
		Short	RK064005	RK065005	RK063005	
22.5" - 25.0" (57.2 cm - 63.5 cm)	L	STD	RK061007	RK062007	RK060007	
		Short	RK064007	RK065007	RK063007	
25.0" - 27.5" (63.5 cm - 69.9 cm)	XL	STD	RK061009	RK062009	RK060009	
		Short	RK064009	RK065009	RK063009	
27.5" - 29.0" (69.9 cm - 73.7 cm)	2XL	STD	RK061011	RK062011	RK060011	
		Short	RK064011	RK065011	RK063011	

### Crossover - Neoprene

Thigh Circumference			Pull-Up		Front Closure	
			w/o Buttress	w/Buttress	w/o Buttress	w/Buttress
15.0" - 17.5" (38.1 cm - 44.5 cm)	XS	STD	RK020001B--	RK021001B--	RK080001B--	RK081001B--
		Short	RK010001B--	RK011001B--	RK070001B--	RK071001B--
17.5" - 20.0" (44.5 cm - 50.8 cm)	S	STD	RK020003B--	RK021003B--	RK080003B--	RK081003B--
		Short	RK010003B--	RK011003B--	RK070003B--	RK071003B--
20.0" - 22.5" (50.8 cm - 57.2 cm)	M	STD	RK020005B--	RK021005B--	RK080005B--	RK081005B--
		Short	RK010005B--	RK011005B--	RK070005B--	RK071005B--
22.5" - 25.0" (57.2 cm - 63.5 cm)	L	STD	RK020007B--	RK021007B--	RK080007B--	RK081007B--
		Short	RK010007B--	RK011007B--	RK070007B--	RK071007B--
25.0" - 27.5" (63.5 cm - 69.9 cm)	XL	STD	RK020009B--	RK021009B--	RK080009B--	RK081009B--
		Short	RK010009B--	RK011009B--	RK070009B--	RK071009B--
27.5" - 29.0" (69.9 cm - 73.7 cm)	2XL	STD	RK020011B--	RK021011B--	RK080011B--	RK081011B--
		Short	RK010011B--	RK011011B--	RK070011B--	RK071011B--



## KNEE SUPPORT

- > Neoprene material provides support, compression, padding, and warmth to the knee
- > Open Patella and Closed Patella versions:
  - > Open Patella: Sewn-in horseshoe buttress
  - > Closed Patella: No buttress

### Indications:

Sore, strained or sprained knee ligaments. Open patella version also indicated for patella tracking support.

Thigh Circumference	Open Patella		Closed Patella	
	Open Patella	Closed Patella	Open Patella	Closed Patella
13.5" - 15.5" (34 cm - 39 cm)	"12" - 15" (30 cm - 38 cm)	XS	SA616001	SA609001
15.5" - 18" (39 cm - 46 cm)	15" - 18" (38 cm - 46 cm)	S	SA616003	SA609003
18" - 20.5" (46 cm - 52 cm)	18" - 21" (46 cm - 53 cm)	M	SA616005	SA609005
20.5" - 23" (52 cm - 58 cm)	21" - 24" (53 cm - 61 cm)	L	SA616007	SA609007
23" - 25.5" (58 cm - 65 cm)	24" - 27" (61 cm - 69 cm)	XL	SA616009	SA609009
25.5" - 28" (65 cm - 71 cm)	27" - 30" (69 cm - 79 cm)	2XL	SA616011	SA609011

Suggested HCFA Code: L1825





## PATELLOFEMORAL STABILIZATION




20.50

REINFORCED PATELLA  
STABILIZER

LATERAL PATELLA  
STABILIZER

TENDON STRAP





## BLED SOE PRINCIPLE:

### ZERO-BASED DESIGN

**B**ledsoe's philosophy of 'zero-based design' looked at Patellofemoral bracing in a whole new light. While most patellofemoral braces function in extension, patellar pain and dislocation are a flexion condition.



RELAXES IN EXTENSION



ACTIVATES IN FLEXION

**The Result:** The first patellofemoral brace to function in flexion.

## 20.50 PATELLOFEMORAL

- > Dynamic control in flexion to control patellar subluxation, dislocation and improve patellofemoral tracking
- > Lightweight Magnesium and short 12" length
- > Also provides additional MCL and LCL support

### Indications:

Patellofemoral Pain Syndrome; patella alta and baja; patellofemoral malalignment, subluxation and dislocation; patellar tendonitis; post operative effusion; jumper's knee; and tight iliotibial band or quadriceps laxity.

- Custom: Fabricated from FitKit 360° measurements, laser scan, or cast mold.
- Off-the-Shelf: Measure 6" above and 6" below mid-patella.

Thigh Circumference	Calf Circumference		Left	Right
13.5" - 16.0" (34.3 cm - 40.6 cm)	11.0" - 12.5" (27.9 cm - 31.7 cm)	XS	PF010101	PF010201
16.0" - 18.75" (40.6 cm - 47.6 cm)	12.5" - 14.0" (31.7 cm - 35.6 cm)	S	PF010103	PF010203
18.75" - 21.5" (47.6 cm - 54.6 cm)	14.0" - 15.5" (35.6 cm - 39.4 cm)	M	PF010105	PF010205
21.5" - 24.25" (54.6 cm - 61.6 cm)	15.5" - 17.0" (39.4 cm - 43.2 cm)	L	PF010107	PF010207
24.25" - 27.0" (61.6 cm - 68.6 cm)	17.0" - 18.5" (43.2 cm - 47.0 cm)	XL	PF010109	PF010209
27.0" - 29.5" (68.6 cm - 74.9 cm)	18.5" - 20.0" (47.0 cm - 50.8 cm)	2XL	PF010111	PF010211

Suggested HCFA Code: L1846 + L2800 (Custom), L1845 + L2795 (OTS)





## REINFORCED PATELLA STABILIZER

- > Neoprene material provides support, padding, compression and warmth
- > Sewn-in tubular horseshoe buttress
- > Adjustable thigh/calf straps provide superior stability and compression

### Indications:

To stabilize and improve patella tracking.

Thigh Circumference		
15.5" - 18" (39 cm - 46 cm)	S	SA614003
18" - 20.5" (46 cm - 52 cm)	M	SA614005
20.5" - 23" (52 cm - 58 cm)	L	SA614007
23" - 25.5" (58 cm - 65 cm)	XL	SA614009
25.5" - 28" (65 cm - 71 cm)	2XL	SA614011

Suggested HCFA Code: L1800 + L2795



## LATERAL PATELLA STABILIZER

- > Neoprene material provides support, padding, compression and warmth
- > Tubular lateral "J" buttress for patellar stabilization
- > Angular strapping for added stability

### Indications:

Patellofemoral disorders including lateral patellar subluxation, dislocations, and lateral patellofemoral malalignment.

Thigh Circumference		Left	Right
13" - 15.5" (33 cm - 39 cm)	XS	SA611101	SA611201
15.5" - 18.5" (39 cm - 47 cm)	S	SA611103	SA611203
18.5" - 21" (47 cm - 53 cm)	M	SA611105	SA611205
21" - 23.5" (53 cm - 60 cm)	L	SA611107	SA611207
23.5" - 26.5" (60 cm - 67 cm)	XL	SA611109	SA611209
26.5" - 29.5" (67 cm - 75 cm)	2XL	SA611111	SA611211

Suggested HCFA Code: L1800 + L2795



## TENDON STRAP

- > Buttress eases the tension of chronic inflammation
- > May also be used on the arm to relieve the pain of 'tennis elbow'

### Indications:


Patella tendonitis and chondromalacia

Knee Circumference		
10" - 11" (25 cm - 28 cm)	XS	SA605001
11" - 12" (28 cm - 30 cm)	S	SA605003
12" - 13" (30 cm - 33 cm)	M	SA605005
13" - 15" (33 cm - 38 cm)	L	SA605007
15" - 17" (38 cm - 46 cm)	XL	SA605009

Suggested HCFA Code: L1825 (Knee), L3701 (Elbow)








2009™  
**Bledsoe**



# POST-OP



EXTENDER PLUS  
MERIT OR  
ORIGINAL  
LEVER LOCK  
CROSSOVER ROM  
THREE PANEL KNEE  
IMMOBILIZER  
SINGLE PANEL KNEE  
IMMOBILIZER  
PHILIPPON HIP  
SIMPLE PELVIC






## BLED SOE PRINCIPLE:

### EASE OF USE

**T**wenty-seven years ago, in the early development of hinged casting, Dr. Poindexter tasked Gary Bledsoe to design a simple removable brace. In his garage, Gary created the industry's first Post-Op Knee Brace. The rest is history...

Today, Bledsoe, a market leader and innovator in Post-Op bracing, applies those same features:

- > Malleable uprights
- > Attached straps to prevent pistoning and improve suspension
- > Breathable one-way stretch foam
- > Adjustable hinge controls flexion and extension from  $-10^{\circ}$  to  $110^{\circ}$



## EXTENDER PLUS

- > Versatile and quick application with push button length adjustment from 22"-34"
- > Accommodates circumference up to 34"
- > Quick release buckles and wrap around design
- > Double padded ankle plates for maximum suspension

### Indications:

Locked or limited motion control of knee during rehabilitation, after operative procedures or injury to the knee ligaments, patellar ligament, articular cartilage, meniscus, and stable or internally fixed fractures of tibial plateau, condyles, proximal tibia or distal femur.

Extender Plus	STD	w/Cuff
Universal	U	EK087000 EK088000

Suggested HCFA Code: L1832

### Shoe Insert

Men's	Women's	Euro	Left	Right	
Up to 4	Up to 5	Up to 34.5	S	EK007101	EK007201
5 - 9	6 - 10	35 - 41	M	EK007105	EK007205





## MERIT OR

- > Optimal suspension strapping configuration
- > Full foam design with four sizes to fit any patient
- > Cool version available for warmer climates

### Indications:

Locked or limited motion control of knee during rehabilitation, after operative procedures or injury to knee ligaments, patellar ligament, articular cartilage, meniscus, and stable or internally fixed fractures of tibial plateau, condyles, proximal tibia or distal femur.

#### Leg Length

16" - 22" (40.6 cm - 55.9 cm)	Short	STD Cool	EK080001 EK083001
22" - 26" (55.9 cm - 66 cm)	Reg	STD Cool	EK080005 EK083005
26" - 30" (66 cm - 76.2 cm)	Long	STD Cool	EK080007 EK083007
30" + (73.7 cm +)	Xlong	STD Cool	EK080009 EK083009

Suggested HCFA Code: L1832



## ORIGINAL

- > 9 sizes for optimal fit
- > Optional posterior and anterior fracture plates
- > Easily converts to KAFO and HKAFO with Bledsoe's Hi-top Boot and/or Simple Pelvic Braces

### Indications:

Locked or limited motion control of knee during rehabilitation, after operative procedures or injury to knee ligaments, patellar ligament, articular cartilage, meniscus, and stable or internally fixed fractures of tibial plateau, condyles, proximal tibia or distal femur. Also indicated for immediate post-injury application for stable, non-displaced or internally fixed fractures of the middle to distal femur and the proximal to middle tibia or fibula when using the optional fracture bracing plates and shoe insert.



#### Shoe Insert

Shoe Size		Euro		Left		Right	
Men's	Women's						
Up to 4	Up to 5	Up to 34.5	S	AM005101-B-	AM005201-B-		
5 - 9	6 - 10	35 - 41	M	AM005105-B-	AM005205-B-		
10 - 14	11 - 15	41.5 - 47.5	L	AM005109-B-	AM005209-B-		

Leg Length	Original Brace	Hinge & Thigh Attachment*		Fracture Plates		
				Left	Right	
17" - 19" (43 cm - 48 cm)	18"	EK020018	EK025018	Femoral Tibial Winged Tibial	AJ005101 AL005101 AL004101	AJ005201 AL005201 AL004201
19" - 21" (48 cm - 53 cm)	20"	EK020020	EK025020	Femoral Tibial Winged Tibial	AJ005101 AL005101 AL004101	AJ005201 AL005201 AL004201
21" - 23" (53 cm - 57 cm)	22"	EK020022	EK025022	Femoral Tibial Winged Tibial	AJ005101 AL005101 AL004101	AJ005201 AL005201 AL004201
23" - 25" (57 cm - 63 cm)	24"	EK020024	EK025024	Femoral Tibial Winged Tibial	AJ005101 AL005101 AL004101	AJ005201 AL005201 AL004201
25" - 27" (63 cm - 69 cm)	26"	EK020026	EK025026	Femoral Tibial Winged Tibial	AJ005105 AL005105 AL004105	AJ005205 AL005205 AL004205
27" - 29" (69 cm - 74 cm)	28"	EK020028	EK025028	Femoral Tibial Winged Tibial	AJ005105 AL005105 AL004105	AJ005205 AL005205 AL004205
29" - 31" (74 cm - 79 cm)	30"	EK020030	EK025030	Femoral Tibial Winged Tibial	AJ005109 AL005109 AL004109	AJ005209 AL005209 AL004209
31" - 33" (79 cm - 84 cm)	32"	EK020032	EK025032	Femoral Tibial Winged Tibial	AJ005109 AL005109 AL004109	AJ005209 AL005209 AL004209
33" - 35" (84 cm - 89 cm)	34"	EK020034	EK025034	Femoral Tibial Winged Tibial	AJ005109 AL005109 AL004109	AJ005209 AL005209 AL004209

Suggested HCFA Code: L1832

\*For use with the Bledsoe Hi-Top Boot



## LEVER LOCK

- > Economy version
- > Easily locks in full extension to immobilize the leg and unlocks to allow motion for rehab, passive exercise, or sitting
- > Ideal for arthroscopic procedures

### Indications:

For application following arthroscopic reconstruction of the ACL or other procedures or injuries which require locking of the knee in full extension with intermediate periods of full range of motion.

Leg Length		
21" - 25" (53.3 cm - 63.5 cm)	Reg	AK070005B
25" - 29" (63.5 cm - 73.7 cm)	Long	AK070007B
29"+ (73.7 cm +)	XLong	AK070009B

Suggested HCFA Code: L1832



## CROSSOVER ROM

- > Cool breathable Spacer Mesh provides support, compression and warmth to the knee
- > Easily adjustable ROM hinges for ideal anatomical placement
- > Front closure design for easy application
- > Interchangeable Sport Hinges convert Crossover ROM to Crossover

### Indications:

Locked or limited motion control during rehabilitation after arthroscopic procedures or injury.

Thigh Circumference		Crossover ROM	Sport Hinges*
15.0" - 17.5" (38.1 cm - 44.5 cm)	XS	RK090001	RK066301
17.5" - 20.0" (44.5 cm - 50.8 cm)	S	RK090003	RK066303
20.0" - 22.5" (50.8 cm - 57.2 cm)	M	RK090005	RK066305
22.5" - 25.0" (57.2 cm - 63.5 cm)	L	RK090007	RK066307
25.0" - 27.5" (63.5 cm - 69.9 cm)	XL	RK090009	RK066309
27.5" - 29.0" (69.9 cm - 73.7 cm)	2XL	RK090011	RK066311

\*Hinges and straps only.

To be used in converting the Crossover ROM to a Crossover.





## THREE PANEL KNEE IMMOBILIZER

- Fits and shapes to legs using medial and lateral stays
- Circumference adjusts up to 30" thigh
- Terry cloth interior for added patient comfort

### Indications:

For motion control of the knee after malalignment injury, post operative care, or during rehabilitation of knee ligaments, cartilage or stable or internally fixed fractures of the tibial plateau, condyles or proximal tibia and distal femur.

### Length

12" (30 cm)	U	SA608012
16" (41 cm)	U	SA608016
20" (51 cm)	U	SA608020
24" (61 cm)	U	SA608024

Suggested HCFA Code: L1830



## SINGLE PANEL KNEE IMMOBILIZER

- Lightweight breathable foam for added patient comfort
- Circumference adjusts up to 30" thigh
- Elastic straps reduce migration

### Indications:

For motion control of the knee after malalignment injury, post operative care, or during rehabilitation of knee ligaments, cartilage or stable or internally fixed fractures of the tibial plateau, condyles or proximal tibia and distal femur.

### Length

12" (30 cm)	U	SA607012
18" (46 cm)	U	SA607018
20" (51 cm)	U	SA607020
22" (56 cm)	U	SA607022
24" (61 cm)	U	SA607024

Suggested HCFA Code: L1830





## PHILIPPON HIP

- > Designed to allow controlled hip flexion in 15° increments from -30° to 105°
- > Can be locked in any position from -30° to 60°
- > Optional KAF Positioning Kit and Post Pad protects patients during and after hip surgery



### Indications:

Post operative application following arthroscopic hip repairs, total hip revisions, or other hip joint surgeries, injuries or problems that can benefit from range-of-motion control.

Patient Height			Left	Right
Under 5' 4" (Under 163 cm)	Short	Brace Pkg	AT090115 AT093115	AT090215 AT093215
4" and Taller (163 cm and Taller)	Reg	Brace Pkg	AT090105 AT093105	AT090205 AT093205

Suggested HCFA Code: L1686



KAF POSITIONING KIT

## SIMPLE PELVIC

- > Designed to allow controlled hip flexion in 15° increments from -30° to 105°
- > Can be locked in 15° increments from -30° to 60°
- > Can be formed for slight abduction
- > Pelvic Attachment easily converts to HKAFO with Original Post-Op Knee and Hi-Top Boot

### Indications:

Post operative application following arthroscopic hip repairs, total hip revisions, or other hip joint surgeries, injuries or problems that can benefit from range-of-motion control.

### Simple Pelvic Series

Leg Length	Iliac Circumference	Thigh Circumference		Simple Pelvic Brace		Pelvic Attachment	
				Left	Right	Left	Right
22" - 26" (56 cm - 66 cm)	38" Max (96.5 cm Max)	26" Max (66 cm Max)	S	AT013101	AT013201	AT010101	AT010201
26" - 30" (66 cm - 76 cm)	45" Max (114 cm Max)	30" Max (76 cm Max)	M	AT013103	AT013203	AT010103	AT010203
30" - 34" (76 cm - 86 cm)	55" Max (140 cm Max)	34" Max (86 cm Max)	L	AT013105	AT013205	AT010105	AT010205

Suggested HCFA Code: L1686






## FOOT & ANKLE



AXIOM ANKLE  
ULTIMATE ANKLE  
SIDEKICK  
WRAPTOR ANKLE  
M-BRACE





## BLED SOE PRINCIPLE:

### ANATOMICAL FIT

When designing the AXIOM Ankle brace, *FIT* was key. It needed to be stable, yet flexible; support the ankle without restriction; and eliminate pressure points that exist in other braces.

Our anatomical design maintains ankle congruency and prevents pistoning.

The AXIOM Ankle “disappears” when the game begins.



## AXIOM ANKLE

- > Low/High ankle sprain protection
- > Prevents inversion and eversion
- > Anatomical design and new hybrid plastic eliminates pressure points

### Indications:

Ideal support for all ankle sprains. From acute grade 1 sprains to chronically unstable ankles. Designed to increase support in inversion and eversion and allow uninhibited dorsi-flexion and plantar-flexion.

Men's	Shoe Size Women's	Euro		Left	Right
up to 7½	up to 9	up to 40	S	AS020103	AS020203
8 - 11½	9½ - 13	40.5 - 44.5	M	AS020105	AS020205
12 +	13½ +	45+	L	AS020107	AS020207

Suggested HCFA Code: L1971





## ULTIMATE ANKLE



- > Attaches to the shoe for total ankle control
- > Adjustable and formable medial and lateral uprights
- > Fully adjustable footplate and uprights allow matching of anatomical pivot points
- > Off-the-shelf and semi-custom for optimal control

### Indications:

Inversion and eversion control of chronically unstable ankles, grade 1 and 2, and acutely sprained ankles where dorsi-flexion and planter-flexion is permissible and inflammation is under control.

● Semi-Custom: See measuring chart (*visit [www.BledsoeBrace.com](http://www.BledsoeBrace.com)*)

● Off-the-Shelf: See measuring chart (*visit [www.BledsoeBrace.com](http://www.BledsoeBrace.com)*)

		Left	Right
RED	02	AS130252	AS230252
ORANGE	03	AS130353	AS230353
YELLOW	04	AS150454	AS250454
GREEN	05	AS150555	AS250555
BLUE	06	AS170656	AS270656
PURPLE	07	AS170757	AS270757

Suggested HCFA Code: L1971

## SIDEKICK

- > Attaches to the shoe for total ankle control
- > Adjustable and formable medial upright
- > Fully adjustable footplate and upright allow matching of anatomical pivot points
- > Single upright for easy application

### Indications:

Inversion and eversion control of chronically unstable ankles, grade 1 and 2, and acutely sprained ankles where dorsi-flexion and planter-flexion is permissible and inflammation is under control.



Medial Ankle Bone Height		Left	Right
2.5" - 3 7/8" (6.5 cm - 9.8 cm)	S	AS002413	AS002423
3 3/8" - 4 7/8" (8.5 cm - 12.5 cm)	L	AS002417	AS002427

Suggested HCFA Code: L1971



## WRAPTOR

- > Quick fit inversion and eversion straps provide figure 8 support
- > Cool breathable mesh liner and tough ballistic nylon exterior for long term wear and patient comfort
- > Also available with Speed Laces for quick application and removal

### Indications:

Chronic and acute ankle sprains, abnormal inversion and eversion control, and as protective sportswear.



#### Ankle/Heel Circumference

10" - 11" (25 cm - 28 cm)	XS	STD Speed	SA702501 SA702001
11" - 12" (28 cm - 30 cm)	S	STD Speed	SA702503 SA702003
12" - 13" (30 cm - 33 cm)	M	STD Speed	SA702505 SA702005
13" - 14" (33 cm - 36 cm)	L	STD Speed	SA702507 SA702007
14" - 15" (36 cm - 38 cm)	XL	STD Speed	SA702509 SA702009
15" - 16" (38 cm - 41cm)	2XL	STD Speed	SA702511 SA702011

Suggested HCFA Code: L1902

## M-BRACE

- > Air pads and air/gel pads available
- > Easy to apply
- > Universal size adjusts to fit most patients

### Indications:

To provide support and reduce swelling of mild sprains.

		M-Brace	20 Pack	50 Pack
U	Air	ES040005	ES042005	ES045005
	Gel	ES050005	ES052005	ES055005
	Combo	.....	ES062005	ES065005

Suggested HCFA Code: L1906





## WALKER BOOTS/ SPLINTS



EVERYDAY BOOTS  
PEDIATRIC BOOT  
LONG TERM BOOTS  
DIABETIC BOOTS  
SPECIALTY BOOTS  
PFS  
SPLINTS  
SILICONE HEEL CUPS






## BLED SOE PRINCIPLE:

### DETAIL DRIVEN

It's the small things that count. Bledsoe evaluates all factors when designing our braces. From the heel height and width, to the shape of the rocker bottom at every phase of the gait cycle, Bledsoe created the safest walking boot.



No pelvic tilt. Safer on inclines.  
Easier on hips and knees.

Smoothest gait in the industry.

## ADJUSTAFIT

### EVERYDAY BOOT

- > Aluminum arms easily adjust to accommodate swelling and post operative dressings
- > Lightweight Aluminum foot shell provides strength without weight
- > Wide version available that still features the proprietary Bledsoe narrow heel

#### Indications:

Ankle sprains and fractures, forefoot fractures, forefoot sprains, and stable distal tibia fractures.



#### AdjustaFit

Men's	Shoe Size Women's	Euro		STD	Pre-Assembled	Mid-Calf	
1 - 2.5	1 - 3.5	29 - 33	XS	w/Air w/o Air w/AHS	AL253001 AL053001 AL053051	AL254001 AL054001 AL054051	AL256001 AL056001 AL057001
2.5 - 5.5	3.5 - 6	33 - 37	S	w/Air w/o Air w/AHS	AL253003 AL053003 AL053053	AL254003 AL054003 AL054053	AL256003 AL056003 AL057003
5.5 - 9	6 - 9.5	37 - 41	M	w/Air w/o Air w/AHS	AL253005 AL053005 AL053055	AL254005 AL054005 AL054055	AL256005 AL056005 AL057005
9 - 12.5	9.5 - 13	41 - 46	L	w/Air w/o Air w/AHS	AL253007 AL053007 AL053057	AL254007 AL054007 AL054057	AL256007 AL056007 AL057007
12.5 - 17	13+	46 - 49+	XL	w/Air w/o Air w/AHS	AL253009 AL053009 AL053059	AL254009 AL054009 AL054059	AL256009 AL056009 AL057009

#### AdjustaFit Wide\*

Men's	Shoe Size Women's	Euro		STD	Pre-Assembled	Mid-Calf	
1 - 2.5	1 - 3.5	29 - 33	XS	w/Air w/o Air	AL253601 AL053601	AL254601 AL054601	AL256601 AL056601
2.5 - 5.5	3.5 - 6	33 - 37	S	w/Air w/o Air	AL253603 AL053603	AL254603 AL054603	AL256603 AL056603
5.5 - 9	6 - 9.5	37 - 41	M	w/Air w/o Air	AL253605 AL053605	AL254605 AL054605	AL256605 AL056605
9 - 12.5	9.5 - 13	41 - 46	L	w/Air w/o Air	AL253607 AL053607	AL254607 AL054607	AL256607 AL056607
12.5 - 17	13+	46 - 49+	XL	w/Air w/o Air	AL253609 AL053609	AL254609 AL054609	AL256609 AL056609

Suggested HCFA Code: L4360 or L4386

\*Wide version pre-shaped 1 inch wider than the standard boot.



## PROGAI

### EVERYDAY BOOT

- > Closed toe bootie
- > Optional ankle-heel pad
- > Optional pneumatic ankle support

#### Indications:

Ankle sprains and fractures, forefoot fractures, forefoot sprains, and stable distal tibia fractures.



Men's	Shoe Size Women's	Euro			STD	Mid-Calf	STD w/ EZ Set	Mid-Calf w/ EZ Set
Kid's 12.5 - Men's 2.5	1 - 3.5	29 - 33	XXS	w/Air w/o Air w/AHS*	AL232000 AL032000 -----	-----	-----	-----
1 - 2.5	1 - 3.5	29 - 33	XS	w/Air w/o Air w/AHS*	AL232001BB- AL032001BB- AL032051BB-	AL252001BB- AL052001BB- AL052051BB-	AL274001BB- AL074001BB- AL074051BB-	AL275001BB- AL075001BB- AL075051BB-
2.5 - 5.5	3.5 - 6	33 - 37	S	w/Air w/o Air w/AHS*	AL232003BB- AL032003BB- AL032053BB-	AL252003BB- AL052003BB- AL052053BB-	AL274003BB- AL074003BB- AL074053BB-	AL275003BB- AL075003BB- AL075053BB-
5.5 - 9	6 - 9.5	37 - 41	M	w/Air w/o Air w/AHS*	AL232005BB- AL032005BB- AL032055BB-	AL252005BB- AL052005BB- AL052055BB-	AL274005BB- AL074005BB- AL074055BB-	AL275005BB- AL075005BB- AL075055BB-
9 - 12.5	9.5 - 13	41 - 46	L	w/Air w/o Air w/AHS*	AL232007BB- AL032007BB- AL032057BB-	AL252007BB- AL052007BB- AL052057BB-	AL274007BB- AL074007BB- AL074057BB-	AL275007BB- AL075007BB- AL075057BB-
12.5 - 17	13+	46 - 49	XL	w/Air w/o Air w/AHS*	AL232009BB- AL032009BB- AL032059BB-	AL252009BB- AL052009BB- AL052059BB-	AL274009BB- AL074009BB- AL074059BB-	AL275009BB- AL075009BB- AL075059BB-

Suggested HCFA Code: L4360 or L4386

\*Ankle Heel Strap (AHS)

## WEE WALKER

### PEDIATRIC BOOT

- > Fits children as young as 1 year
- > Heel height and width and rocker bottom designed especially for small children
- > Forefoot straps immobilize the ankle

#### Indications:

Ankle sprains and fractures, forefoot fractures, forefoot sprains, and stable distal tibia fractures.

Age in Years	Foot Length	Leg Length		Wee Walker
1 - 2	up to 6" (15 cm)	up to 9" (23 cm)	S-M	AL185003BB-
2 1/2 - 6	6" - 7.5" (15 cm - 19 cm)	up to 11" (28 cm)	L	AL185007BB-
6 1/2 - 9 1/2	7.5" - 8.5" (19 cm - 22 cm)	up to 14" (35.5 cm)	XL	AL185009BB-

Suggested HCFA Code: L4386





## LONG TERM BOOT



- > Our most durable fracture boots
- > Dual-layer foam for longer wear
- > Optional EZ Set Hinge, and Patella Tendon Bearing and Anterior Fracture Plates available
- > Three heights to meet your specific needs
  - > Hi-Top: 2 inches taller than our standard boot height for maximum control
  - > Lo-Top: Standard boot height
  - > Mid-Calf: Short boot height

### Indications:

Ankle sprains and fractures, forefoot fractures, forefoot sprains, and stable distal tibia fractures. Indications include middle to distal fibula or tibia fractures when fracture plates are used.

### Hi-Top Boot

Men's	Shoe Size		Euro		STD		EZ Set	
	Women's				w/Cuff	w/ Bootie	w/Cuff	w/ Bootie
1 - 2.5	1 - 3.5	29 - 33	XS	w/Air w/o Air	AL236001BB- AL036001BB-	AL236021BB- AL036021BB-	AL272001BB- AL072001BB-	AL272021BB- AL072021BB-
2.5 - 5.5	3.5 - 6	33 - 37	S	w/Air w/o Air	AL236003BB- AL036003BB-	AL236023BB- AL036023BB-	AL272003BB- AL072003BB-	AL272023BB- AL072023BB-
5.5 - 9	6 - 9.5	37 - 41	M	w/Air w/o Air	AL236005BB- AL036005BB-	AL236025BB- AL036025BB-	AL272005BB- AL072005BB-	AL272025BB- AL072025BB-
9 - 12.5	9.5 - 13	41 - 46	L	w/Air w/o Air	AL236007BB- AL036007BB-	AL236027BB- AL036027BB-	AL272007BB- AL072007BB-	AL272027BB- AL072027BB-
12.5 - 17	13+	46 - 49+	XL	w/Air w/o Air	AL236009BB- AL036009BB-	AL236029BB- AL036029BB-	AL272009BB- AL072009BB-	AL272029BB- AL072029BB-

Suggested HCFA Code: L4360 or L4386

### Hi-Top PTB

Men's	Shoe Size		Euro		PTB w/ Cuff		PTB w/ Bootie	
	Women's				Left	Right	Left	Right
1 - 2.5	1 - 3.5	29 - 33	XS	w/Air w/o Air	AL235101BB- AL035101BB-	AL235201BB- AL035201BB-	AL235121BB- AL035121BB-	AL235221BB- AL035221BB-
2.5 - 5.5	3.5 - 6	33 - 37	S	w/Air w/o Air	AL235103BB- AL035103BB-	AL235203BB- AL035203BB-	AL235123BB- AL035123BB-	AL235223BB- AL035223BB-
5.5 - 9	6 - 9.5	37 - 41	M	w/Air w/o Air	AL235105BB- AL035105BB-	AL235205BB- AL035205BB-	AL235125BB- AL035125BB-	AL235225BB- AL035225BB-
9 - 12.5	9.5 - 13	41 - 46	L	w/Air w/o Air	AL235107BB- AL035107BB-	AL235207BB- AL035207BB-	AL235127BB- AL035127BB-	AL235227BB- AL035227BB-
12.5 - 17	13+	46 - 49+	XL	w/Air w/o Air	AL235109BB- AL035109BB-	AL235209BB- AL035209BB-	AL235129BB- AL035129BB-	AL235229BB- AL035229BB-

Suggested HCFA Code: L4360 or L4386

### Lo-Top

Men's	Shoe Size		Euro		STD		EZ Set	
	Women's				w/ Cuff	w/ Bootie	w/ Cuff	w/ Bootie
1 - 2.5	1 - 3.5	29 - 33	XS	w/Air w/o Air	AL285001BB- AL085001BB-	AL285021BB- AL085021BB-	AL270001BB- AL070001BB-	AL270021BB- AL070021BB-
2.5 - 5.5	3.5 - 6	33 - 37	S	w/Air w/o Air	AL285003BB- AL085003BB-	AL285023BB- AL085023BB-	AL270003BB- AL070003BB-	AL270023BB- AL070023BB-
5.5 - 9	6 - 9.5	37 - 41	M	w/Air w/o Air	AL285005BB- AL085005BB-	AL285025BB- AL085025BB-	AL270005BB- AL070005BB-	AL270025BB- AL070025BB-
9 - 12.5	9.5 - 13	41 - 46	L	w/Air w/o Air	AL285007BB- AL085007BB-	AL285027BB- AL085027BB-	AL270007BB- AL070007BB-	AL270027BB- AL070027BB-
12.5 - 17	13+	46 - 49+	XL	w/Air w/o Air	AL285009BB- AL085009BB-	AL285029BB- AL085029BB-	AL270009BB- AL070009BB-	AL270029BB- AL070029BB-

Suggested HCFA Code: L4360 or L4386

### Mid-Calf Boot

Men's	Shoe Size		Euro		STD		EZ Set	
	Women's				w/Air w/o Air	w/Air w/o Air	w/Air w/o Air	w/Air w/o Air
1 - 2.5	1 - 3.5	29 - 33	XS	w/Air w/o Air	AL255001BB- AL055001BB-	AL271001BB- AL071001BB-		
2.5 - 5.5	3.5 - 6	33 - 37	S	w/Air w/o Air	AL255003BB- AL055003BB-	AL271003BB- AL071003BB-		
5.5 - 9	6 - 9.5	37 - 41	M	w/Air w/o Air	AL255005BB- AL055005BB-	AL271005BB- AL071005BB-		
9 - 12.5	9.5 - 13	41 - 46	L	w/Air w/o Air	AL255007BB- AL055007BB-	AL271007BB- AL071007BB-		
12.5 - 17	13+	46 - 49+	XL	w/Air w/o Air	AL255009BB- AL055009BB-	AL271009BB- AL071009BB-		

Suggested HCFA Code: L4360 or L4386




Mid-Calf



## BLEDSOE PRINCIPLE:


OUTPERFORM

 ur Diabetic Boots not only protect the foot they reduce peak pressure--30% to be exact!



A combination of proprietary Auto-mold foam and dual density mid-sole promotes perimeter unloading and prevents movement that can breakdown the skin.

That is why the Bledsoe Conformer is recognized as the best performing Diabetic Boot.



## CONFORMER

DIABETIC BOOT

- Outperforms Total Contact Cast in reducing pressure—30% less peak pressure
- Proprietary Auto-mold™ foam contours to the foot and prevents movement that can breakdown the skin
- Dual density mid-sole foam promotes perimeter loading

### Indications:

Ulcerative/pre-ulcerative conditions on the bottom of the foot for sensate and non-sensate patients.

Men's	Shoe Size Women's	Euro	Size		Left	Right
4 - 5.5	5 - 6.5	35 - 37	04	w/Air	AL241104	AL241204
				w/o Air	AL041104	AL041204
6 - 7.5	7 - 8.5	37.5 - 39.5	06	w/Air	AL241106	AL241206
				w/o Air	AL041106	AL041206
8 - 9.5	9 - 10.5	40 - 42	08	w/Air	AL241108	AL241208
				w/o Air	AL041108	AL041208
10 - 11.5	11 - 12.5	42.5 - 44.5	10	w/Air	AL241110	AL241210
				w/o Air	AL041110	AL041210
12 - 13.5	13 - 14.5	45 - 47	12	w/Air	AL241112	AL241212
				w/o Air	AL041112	AL041212
14 - 15.5	15 - 16.5	47.5 - 49.5	14	w/Air	AL241114	AL241214
				w/o Air	AL041114	AL041214

Suggested HCFA Code: L4360 or L4386





# CHARCOT CONFORMER

## DIABETIC BOOT

- > 2 inches taller than standard Conformer
- > Pneumatic Patella Tendon Bearing adjustment for comfort and optimal offloading

### Indications:

Ulcerative/pre-ulcerative conditions on the bottom of the foot for sensate and non-sensate patients when partial patellar tendon bearing is desired.

Men's	Shoe Size		Ankle Size	Ankle		
	Women's	Euro		Air	Left	Right
4 - 5.5	5 - 6.5	35 - 37	04	w/Air w/o Air	AL242104 AL042104	AL242204 AL042204
6 - 7.5	7 - 8.5	37.5 - 39.5	06	w/Air w/o Air	AL242106 AL042106	AL242206 AL042206
8 - 9.5	9 - 10.5	40 - 42	08	w/Air w/o Air	AL242108 AL042108	AL242208 AL042208
10 - 11.5	11 - 12.5	42.5 - 44.5	10	w/Air w/o Air	AL242110 AL042110	AL242210 AL042210
12 - 13.5	13 - 14.5	45 - 47	12	w/Air w/o Air	AL242112 AL042112	AL242212 AL042212
14 - 15.5	15 - 16.5	47.5 - 49.5	14	w/Air w/o Air	AL242114 AL042114	AL242214 AL042214

Suggested HCFA Code: L4360





## BLED SOE PRINCIPLE:

SPECIALIZE

**B**oots are one of our specialties. We understand that not every patient's needs are the same. That's why we make a variety of boots tailored for specific applications. All still have that attention to detail that Bledsoe is known for when it comes to boots.



Bunion



Platform



Achilles



Nite N' Day

## BUNION

SPECIALTY BOOT

- > Pad-and-band maintains proper alignment of the Great Toe
- > Soft Auto-mold™ sole cushions the foot while preventing unnecessary movement and pressure on the toe
- > Pneumatic Forefoot pad
- > Patients are able to walk sooner without crutches

### Indications:

Post-Operative application for Hallux Abducto Valgus (HAV) surgery.

Men's	Shoe Size Women's	Euro	Ankle Air	
2.5 - 5.5	3.5 - 6	33 - 37	S	w/Air AL258003 w/o Air AL058003
5.5 - 9	6 - 9.5	37 - 41	M	w/Air AL258005 w/o Air AL058005
9 - 12.5	9.5 - 13	41 - 46	L	w/Air AL258007 w/o Air AL058007

Suggested HCFA Code: L4360





## PLATFORM

### SPECIALTY BOOT

- > Flat innersole helps prevent dorsi-flexion for post-operative tarsal, metatarsal and phalange procedures
- > May also be used for metatarsal, forefoot or calcaneus fractures and forefoot sprains

#### Indications:

Post-Op for tarsal, metatarsal, and phalange procedures

Men's	Shoe Size Women's	Euro			STD	Mid-Calf
1 - 2.5	1 - 3.5	29 - 33	XS	w/Air w/o Air	AL243001BB- AL043001BB-	AL243101 AL043101
2.5 - 5.5	3.5 - 6	33 - 37	S	w/Air w/o Air	AL243003BB- AL043003BB-	AL243103 AL043103
5.5 - 9	6 - 9.5	37 - 41	M	w/Air w/o Air	AL243005BB- AL043005BB-	AL243105 AL043105
9 - 12.5	9.5 - 13	41 - 46	L	w/Air w/o Air	AL243007BB- AL043007BB-	AL243107 AL043107
12.5 - 17	13+	46 - 49+	XL	w/Air w/o Air	AL243009BB- AL043009BB-	AL243109 AL043109

Suggested HCFA Code: L4360 or L4386



## ACHILLES

### SPECIALTY BOOT

- > Places the foot in plantar-flexed position, up to 40° using four 10° wedges in the heel—as the condition improves, wedges are removed
- > Gently stretches the Achilles Tendon during rehab without casting
- > Wedges allow the patient to avoid pressure applied to the Achilles Tendon

#### Indications:

For application after injury or operative procedure of the Achilles Tendon.

Men's	Shoe Size Women's	Euro			STD	PTB	
						Left	Right
2.5 - 5.5	3.5 - 6	33 - 37	S	w/Air w/o Air	AL236013BB- AL036013BB-	AL235113BB- AL035113BB-	AL235213BB- AL035213BB-
5.5 - 9	6 - 9.5	37 - 41	M	w/Air w/o Air	AL236015BB- AL036015BB-	AL235115BB- AL035115BB-	AL235215BB- AL035215BB-
9 - 12.5	9.5 - 13	41 - 46	L	w/Air w/o Air	AL236017BB- AL036017BB-	AL235117BB- AL035117BB-	AL235217BB- AL035217BB-
12.5 - 17	13+	46 - 49+	XL	w/Air w/o Air	AL236019BB- AL036019BB-	AL235119BB- AL035119BB-	AL235219BB- AL035219BB-

Suggested HCFA Code: L4360 or L4386





## NITE N' DAY (PLANTAR FASCIITIS BOOT)

### SPECIALTY BOOT

- > A combination of wedges fit beneath the forefoot and toe keeping the foot in a dorsi-flexed position—as the condition improves, wedges are removed
- > Prevents tears in fibrous tissue that occur after a night of rest—breaking the night time healing/morning pain cycle
- > Includes Night Cover

#### Indications:

Plantar Fasciitis.

Men's	Shoe Size Women's	Euro			
1 - 2.5	1 - 3.5	29 - 33	XS	w/Air w/o Air	AL243011 AL043011
2.5 - 5.5	3.5 - 6	33 - 37	S	w/Air w/o Air	AL243013 AL043013
5.5 - 9	6 - 9.5	37 - 41	M	w/Air w/o Air	AL243015 AL043015
9 - 12.5	9.5 - 13	41 - 46	L	w/Air w/o Air	AL243017 AL043017
12.5 - 17	13+	46 - 49+	XL	w/Air w/o Air	AL243019 AL043019

Suggested HCFA Code: L4360 or L4386



## PFS (PLANTAR FASCIITIS STRAP)

- > Reduces Plantar Fasciitis pain through continuous elastic tension and pressure
- > Includes patient exercise program to control pain and strengthen the muscles of the foot, arch and lower leg

#### Indications:

Plantar Fasciitis.



Men's	Shoe Size Women's	Euro			Left	Right
0 - 8	0 - 8.5	29 - 40	S	STD 10 Pk	ES007103 ES008103	ES007203 ES008203
8.5 - 11	9 - 11.5	40.5 - 43.5	M	STD 10 Pk	ES007105 ES008105	ES007205 ES008205
11.5 +	12 +	44 +	L	STD 10 Pk	ES007107 ES008107	ES007207 ES008207

Suggested HCFA Code: L1902



## PADDED PF SPLINT

### SPLINT

- > Padded footplate keeps the foot in a dorsi-flexed position
- > Prevents tears in fibrous tissue that occur after a night of rest, breaking the night time healing/morning pain cycle
- > Soft foam for nighttime wear

#### Indications:

Plantar Fasciitis, Achilles Tendonitis, plantar flexion contractures, and heel spurs.

Men's	Shoe Size		Euro	S	SA902503
	Women's				
5.5 - 6.5	5 - 6	35 - 36.5	S	SA902503	
7 - 8	6.5 - 7.5	37 - 38.5	M	SA902505	
8.5 - 9.5	8 - 9	39 - 40.5	L	SA902507	
10 - 11	9.5 - 10.5	41 - 42.5	XL	SA902509	
11.5 +	11 +	43 +	2XL	SA902511	

Suggested HCFA Code: L4396



## NIGHT SPLINT

### SPLINT

- > Angled polyethylene shell keeps the foot in a dorsi-flexed position while straps tension for support
- > Prevents tears in fibrous tissue that occur after a night of rest, breaking the night time healing/morning pain cycle
- > Washable liner

#### Indications:

Plantar Fasciitis, Achilles Tendonitis, plantar flexion contractures, and heel spurs.

Men's	Shoe Size		Euro	S	SA902003
	Women's				
0 - 7	0 - 7.5	0 - 38.5	S	SA902003	
7.5 - 10	8 - 10.5	39 - 42	M	SA902005	
10.5 +	11 +	42.5 +	L	SA902007	

Suggested HCFA Code: L4396





## SILICONE HEEL CUPS


- > Dual-density medical-grade silicone to absorb heel pressure
- > Soft silicone center cushions sensitive areas
- > Use in dress or athletic shoes

### Indications:

For treatment of heel spurs, leg pain, general foot fatigue and for treatment of Plantar Fasciitis or Achilles Tendonitis.

Shoe Size		Euro	Pair
Men's	Women's		
3.5 - 6.5	4 - 7	34 - 38	S SA903003
7 - 10.5	7.5 - 11	38.5 - 43.5	M SA903005
11+	11.5+	44+	L SA903007

Suggested HCFA Code: L3485



## POST-OP SHOE

### SPLINT

- > Protects foot after surgery
- > Trimable straps for easy application
- > Large opening to accommodate dressings and swelling

### Indications:

For use following surgery or trauma to the foot.

Shoe Size		Men's	Women's	Pair
Men's	Women's			
6 - 8	4 - 6	S	Men's Women's	SA900003 SA900503
8 - 10	6 - 8	M	Men's Women's	SA900005 SA900505
10 - 12	8 - 10	L	Men's Women's	SA900007 SA900507

Suggested HCFA Code: L3260



## CAST SHOE

### SPLINT

- > Protects cast from dirt and wear
- > No-skid sole
- > Universal right or left

### Indications:


For lower extremity casted fractures.

Length	Width		
7.25" (18 cm)	3.5" (9 cm)	XS	SA901001
9.2" (24 cm)	4.2" (11 cm)	S	SA901003
11" (28 cm)	4.75" (12 cm)	M	SA901005
12.5" (32 cm)	5" (13 cm)	L	SA901007
13.5" (34 cm)	6" (15 cm)	XL	SA901009

Suggested HCFA Code: L3260








2009™  
**Bledsoe**



## UPPER EXTREMITY



ARC / ARC - XR  
ARC WITH PILLOW  
CLINIC ARM SLING  
CLINIC SHOULDER  
IMMOBILIZER  
UNIVERSAL  
ARM SLING  
SLING AND SWATHE  
IMMOBILIZER  
CLAVICLE SUPPORT





## BLED SOE PRINCIPLE:

A NEW PARADIGM

**B**ledsoe challenged its engineers to come up with something so innovative it would change the way the industry viewed shoulder bracing. The result, the first advancement in shoulder bracing in 12 years! A no-pillow design that reliably HOLDS position while reducing SKU's from 16 to 2.



INNOVATIVE ARC FRAME  
POSITIONING SYSTEM

## ARC & ARC-XR

- > ARC ranges from 20° internal rotation to 15° external rotation
- > ARC-XR allows additional rotation from 15° to 70° external rotation
- > Aluminum waistband molds to the patient's torso and holds its position with no awkward pillow
- > One brace fits all patients—right or left, football player or gymnast

### Indications:

For non-surgical or post operative support and immobilization of the injured shoulder where glenohumeral positioning is key.

	ARC	ARC-XR
U STD	AE050100	AE050200
Cool	AE055100	AE055200

**Bledsoe**  
**ARC**  
SHOULDER BRACE





## ARC W/ PILLOW

**Bledsoe**  
**ARC**  
SHOULDER BRACE


- > One brace fits all patients—right or left, football player or gymnast
- > Traditional pillow design

### Indications:

For non-surgical or post operative support and immobilization of the injured shoulder where glenohumeral positioning is key.

U	STD	AE050300
U	Cool	AE055300

Suggested HCFA Code: L3670



## CLINIC ARM SLING

- > Over-sized envelope accommodates large casts and dressings
- > Customized fitting and adjustment feature on front and back
- > Padded strap with easy touch closure

### Indications:

Immobilization and support of mild sprains and strains of the shoulder or fractures of the arm.

Envelope Length		Clinic Arm Sling
11" (28 cm)	XS	SA200001
13.5" (34 cm)	S	SA200003
15.5" (39 cm)	M	SA200005
17.5" (44 cm)	L	SA200007
20" (51 cm)	XL	SA200009

Suggested HCFA Code: A4565



## CLINIC SHOULDER IMMOBILIZER

- > Immobilizing strap secures shoulder in place
- > Over-sized envelope accommodates large casts and dressings
- > Customized fitting and adjustment feature on front and back
- > Padded strap with easy touch closure

### Indications:

Immobilization and support of mild sprains and strains of the shoulder or fractures of the arm.

Envelope Length		
11" (28 cm)	XS	SA200501
13.5" (34 cm)	S	SA200503
15.5" (39 cm)	M	SA200505
17.5" (44 cm)	L	SA200507
20" (51 cm)	XL	SA200509

Suggested HCFA Code: L3670





## UNIVERSAL ARM SLING

- > No trimming necessary, simply fold to proper length to accommodate most patients
- > Large envelope accommodates casts and dressings
- > Easily adjustable slide buckle closure

### Indications:

Immobilization and support of mild sprains and strains of the shoulder or fractures of the arm.



Universal Arm Sling

U	SA201000
---	----------

Suggested HCFA Code: A4565

## SLING AND SWATHE IMMOBILIZER

- > Lightweight, soft foam sling supports the arm
- > Swathe immobilizes the arm and shoulder
- > Adjustable slide buckle straps

### Indications:

Immobilization and support of mild sprains and strains and fractures. Swathe firmly secures arm against body or fixates shoulder.



Universal Arm Sling

U	SA204000
---	----------

Suggested HCFA Code: L3670

## CLAVICLE SUPPORT

- > Urethane foam contoured for fit and comfort
- > Easy touch closures are fully adjustable
- > Easily adjust tension with rear buckles

### Indications:

Clavicular fractures and posture correction.

### Clavicle Support

Chest Circumference		Clavicle Support
20" - 24" (51 cm - 61 cm)	XS	SA301001
24" - 30" (61 cm - 76 cm)	S	SA301003
30" - 36" (76 cm - 91 cm)	M	SA301005
36" - 42" (91 cm - 107 cm)	L	SA301007
42" - 48" (107 cm - 122 cm)	XL	SA301009

Suggested HCFA Code: L3650





## ELBOW, WRIST AND HAND



TELESCOPING ELBOW

TENDON STRAP


WRIST LACER

WRIST LACER WITH  
THUMB SPICA

WRIST BRACE

WRIST BRACE WITH  
THUMB SPICA

TRIGGER AND KEEPER  
THUMB SUPPORT





## TELESCOPING ELBOW

- > Push-button telescoping and malleable arm plates for perfect fit
- > Locks the elbow in any position from  $-10^{\circ}$  to  $110^{\circ}$  in  $10^{\circ}$  increments
- > Easy to apply—no tools required

### Indications:

For immediate application following injury or operative procedures to the elbow, or for stable fractures of the distal humerus or proximal to middle radius or ulna.

		Left	Right
Telescoping Elbow Brace	U	AE023100	AE023200
Supination/Pronation Hand Attachment	U	AE012123B--	AE012223B--
Neutral Hand Attachment	S	AE012111B--	AE012211B--
	L	AE012113B--	AE012213B--

Suggested HCFA Code: L3760



## TENDON STRAP

- > Buttress eases the pain of 'tennis elbow'
- > May also be used on the leg to relieve tension of chronic patellar inflammation

### Indications:

Tennis Elbow, Carpal Tunnel Syndrome and other forms of forearm and wrist pain.

Arm Circumference		
10" - 11" (25 cm - 28 cm)	XS	SA605001
11" - 12" (28 cm - 30 cm)	S	SA605003
12" - 13" (30 cm - 33 cm)	M	SA605005
13" - 15" (33 cm - 38 cm)	L	SA605007
15" - 17" (38 cm - 46 cm)	XL	SA605009

Suggested HCFA Code: L1825 (Knee), L3701 (Elbow)





## WRIST LACER

- > Malleable palmer stays provide customized angulations
- > Thumb cut-out improves patient manipulation, movement and compliance
- > Durable suede is perforated for cool comfort

### Indications:

Immobilization and support for Carpel Tunnel Syndrome, rehab and post-cast support.

Wrist Circumference		Left	Right
Up to 5" (14 cm)	XS	8" SA800101 10" SA800151	SA800201 SA800251
5" - 6" (14 cm - 16 cm)	S	8" SA800103 10" SA800153	SA800203 SA800253
6" - 7" (16 cm - 18 cm)	M	8" SA800105 10" SA800155	SA800205 SA800255
7" - 8" (18 cm - 20 cm)	L	8" SA800107 10" SA800157	SA800207 SA800257
8" - 9" (20 cm - 22 cm)	XL	8" SA800109 10" SA800159	SA800209 SA800259

Suggested HCFA Code: L3908



## WRIST LACER WITH THUMB SPICA

- > Malleable palmer stays provide customized positioning
- > Extensor surface of thumb
- > Durable suede is perforated for cool comfort

### Indications:

For support and symptomatic relief of Gamekeeper's Thumb, deQuervain's Syndrome, scaphoid injury, Carpel Tunnel, and sprains and strains.

Wrist Circumference		Left	Right
Up to 5" (14 cm)	XS	SA805101	SA805201
5" - 6" (14 cm - 16 cm)	S	SA805103	SA805203
6" - 7" (16 cm - 18 cm)	M	SA805105	SA805205
7" - 8" (18 cm - 20 cm)	L	SA805107	SA805207
8" - 9" (20 cm - 22 cm)	XL	SA805109	SA805209

Suggested HCFA Code: L3807



## WRIST BRACE

- > Aluminum stay pre-formed for anatomically correct fit
- > Comfortable, lightweight foam laminate

### Indications:

Reduces the symptoms of Carpal Tunnel Syndrome and effectively supports sprains and strains.

Wrist Circumference		Left	Right
Up to 5" (14 cm)	XS	8" SA802101 10" SA802151	SA802201 SA802251
5" - 6.5" (14 cm - 17 cm)	S	8" SA802103 10" SA802153	SA802203 SA802253
6.5" - 7.5" (17 cm - 19 cm)	M	8" SA802105 10" SA802155	SA802205 SA802255
7.5" - 8.5" (19 cm - 22 cm)	L	8" SA802107 10" SA802157	SA802207 SA802257
8.5" + (22 cm +)	XL	8" SA802109 10" SA802159	SA802209 SA802259

Suggested HCFA Code: L3908



## WRIST BRACE WITH THUMB SPICA

- > Dual Aluminum stays for customized positioning
- > Palmer surface of wrist
- > Extensor surface of thumb

### Indications:

For support and symptomatic relief of Gamekeeper's Thumb, deQuervain's Syndrome, scaphoid injury, Carpel Tunnel, and sprains and strains.

Wrist Circumference		Left	Right
Up to 5" (14 cm)	XS	SA803101	SA803201
5" - 6" (14 cm - 16 cm)	S	SA803103	SA803203
6" - 7" (16 cm - 18 cm)	M	SA803105	SA803205
7" - 8" (18 cm - 20 cm)	L	SA803107	SA803207
8" - 9" (20 cm - 22 cm)	XL	SA803109	SA803209

Suggested HCFA Code: L3807





## TRIGGER AND KEEPER THUMB SUPPORT

---

- > Wrap around one-hand application
- > Flexible support for customized positioning
- > Fits right and left hand

### Indications:

For immobilization, support and relief of Gamekeeper's Thumb, Basal Joint Arthritis, and Thumb Tendonitis.

Wrist Circumference		T & K
Up to 6.5" (17 cm)	S	SA801003
6.5" - 7.5" (17 cm - 19 cm)	M	SA801005
7.5" + (19 cm +)	L	SA801007

Suggested HCFA Code: L3923



2009™  
© Bledsoe



# SPINE




VOLARE W/CHAIRBACK

VOLARE LS-XT

VOLARE LS

VOLARE LITE



**Bledsoe**™



## BLED SOE PRINCIPLE:

### FLEXIBILITY

**W**hether sitting, walking, or lifting, the Volare instantly adjusts for activity with a quick-pull zip line. Independent guides also adjust for girth and shape and allow coning for wide variations in hip-to-waist measurements. As a result, patients are supported throughout the day and throughout their recovery.



## VOLARE WITH CHAIRBACK

- > Anterior and Posterior panels unload spinal pressure for pain relief
- > Zip Line provides 5:1 mechanical advantage for easy one-hand adjustment
- > Posterior chairback panel extends from sacrococcygeal junction to T-9 vertebra

### Indications:

Mechanical back pain, Posterior Lateral Fusion, Lumbar Spinal Laminectomy, multiple level decompression, Spondylolysis, Spondylolisthesis, and Spinal Stenosis.

Waist				
25" - 30" (63.5 cm - 76 cm)	S	8"	10"	AB027003 AB028003
30" - 35" (76 cm - 89 cm)	M	8"	10"	AB027005 AB028005
35" - 40" (89 cm - 101.5 cm)	L	8"	10"	AB027007 AB028007
40" - 45" (101.5 cm - 114 cm)	XL	8"	10"	AB027009 AB028009
45" - 50" (114 cm - 127 cm)	2XL	8"	10"	AB027011 AB028011
50" - 55" (127 cm - 140 cm)	3XL	8"	10"	AB027013 AB028013
56" - 68" (140 cm - 173 cm)	4XL	8"	10"	..... AB028015



PDAC Approved: L0637



## VOLARE LS-XT

- > Anterior and Posterior panels unload spinal pressure for pain relief
- > Zip Line provides 5:1 mechanical advantage for easy one-hand adjustment
- > Posterior panel extends from sacrococcygeal junction to T-9 vertebra
- > Anterior and posterior control

### Indications:

Thoracolumbar injury, revision surgery, multi-level fusion, mechanical back pain, disc herniation, disc degeneration, and Lumbar Discectomy.

### Waist

25" - 30" (63.5 cm - 76 cm)	S	8" 10"	AB022003 AB023003
30" - 35" (76 cm - 89 cm)	M	8" 10"	AB022005 AB023005
35" - 40" (89 cm - 101.5 cm)	L	8" 10"	AB022007 AB023007
40" - 45" (101.5 cm - 114 cm)	XL	8" 10"	AB022009 AB023009
45" - 50" (114 cm - 127 cm)	2XL	8" 10"	AB022011 AB023011
50" - 55" (127 cm - 140 cm)	3XL	8" 10"	AB022013 AB023013
56" - 68" (140 cm - 173 cm)	4XL	8" 10"	----- AB023015

PDAC Approved: L0631



## VOLARE LS

- > Anterior and Posterior panels unload spinal pressure for pain relief
- > Zip Line provides 5:1 mechanical advantage for easy one-hand adjustment
- > Posterior panel extends from L-1 to below S-1 vertebra
- > Anterior and posterior control

### Indications:

Multiple level decompression, Lumbar Laminectomy, Posterior Lateral Fusion, and mechanical back pain.

### Waist

25" - 30" (63.5 cm - 76 cm)	S	8" 10"	AB017003 AB018003
30" - 35" (76 cm - 89 cm)	M	8" 10"	AB017005 AB018005
35" - 40" (89 cm - 101.5 cm)	L	8" 10"	AB017007 AB018007
40" - 45" (101.5 cm - 114 cm)	XL	8" 10"	AB017009 AB018009
45" - 50" (114 cm - 127 cm)	2XL	8" 10"	AB017011 AB018011
50" - 55" (127 cm - 140 cm)	3XL	8" 10"	AB017013 AB018013

PDAC Approved: L0627





## VOLARE LITE

- > Relieves day-to-day back pain from activity, sports or simple desk fatigue
- > Zip Line provides 3:1 mechanical advantage for easy one-hand adjustment

### Indications:

Lower back pain, spinal instability, industrial back support, back sprains and strains, and mechanical back pain.

Waist			
25" - 30"	(63.5 cm - 76 cm)	S	AB012003
30" - 35"	(76 cm - 89 cm)	M	AB012005
35" - 40"	(89 cm - 101.5 cm)	L	AB012007
40" - 45"	(101.5 cm - 114 cm)	XL	AB012009
45" - 50"	(114 cm - 127 cm)	2XL	AB012011

PDAC Approved: L0626



## MODULAR SYSTEM



Low Anterior Panel

AB007002



xt Posterior Panel

AB007006



Mid Anterior Posterior Panel

AB007004

### Lordotic Pads

Correction	
7°	AB005001
10°	AB005002
12°	AB005003
15°	AB005004
10, 12°	AB005005
7, 10, 12, 15°	AB005006



### Chairback Frame

Waist		
25" - 30"	S	AB008005
(63.5 cm - 76 cm)		
30" - 35"	M	AB008005
(76 cm - 89 cm)		
35" - 40"	L	AB008007
(89 cm - 101.5 cm)		
40" - 45"	XL	AB008009
(101.5 cm - 114 cm)		
45" - 50"	2XL	AB008011
(114 cm - 127 cm)		
50" - 55"	3XL	AB008011
(127 cm - 140 cm)		
56" - 68"	4XL	AB008011
(140 cm - 173 cm)		



### Volare Corset

Waist				
25" - 30"	S	8"	AB013003	
(63.5 cm - 76 cm)		10"	AB014003	
30" - 35"	M	8"	AB013005	
(76 cm - 89 cm)		10"	AB014005	
35" - 40"	L	8"	AB013007	
(89 cm - 101.5 cm)		10"	AB014007	
40" - 45"	XL	8"	AB013009	
(101.5 cm - 114 cm)		10"	AB014009	
45" - 50"	2XL	8"	AB013011	
(114 cm - 127 cm)		10"	AB014011	
50" - 55"	3XL	8"	AB013013	
(127 cm - 140 cm)		10"	AB014013	





## COLD THERAPY



COLD CONTROL

UNIVERSAL PAD

SHOULDER PAD

HAND/WRIST/FINGERS PAD

ANKLE PAD

5X10 RECTANGLE PAD

10X14 RECTANGLE PAD

SHOULDER WRAP





## BLED SOE PRINCIPLE:

SAFETY FIRST

**B**ledsoe didn't invent cold therapy—we perfected it!

Patient safety is critical; even when dealing with the safest anti-inflammatory medication—ice. The Cold Control eliminates the risk of nerve damage or frostbite.

Patients are protected while sleeping or from prolonged therapy.



## COLD CONTROL

- > Eliminates the risk of nerve damage or frostbite
- > Patient Temperature Control
- > Greater ice capacity holding as little as 6 quarts or as much as 14 quarts for longer continuous therapy

### Indications:

For immediate application following injury or operative procedures to reduce pain and swelling and aid in healing.

			Single	6 Pack
Unit			SV005100	SV105100
w/ Universal Pad	Regular Hose	Non-Sterile	SV007002	SV107002
w/ Universal Pad	Extended Hose	Non-Sterile	SV007102	-----
w/ Shoulder Pad	Extended Hose	Non-Sterile	SV007106	-----






## COLD CONTROL PADS

### UNIVERSAL

- > Universal shoulder, elbow, knee, and ankle
- > Regular hose lengths and extended hose lengths available

			Single	20 Pack
Regular Hose	U	Non-Sterile	SV003002	SV103002
Regular Hose	U	Sterile	SV013002S	SV113002S
Extended Hose	U	Non-Sterile	SV003102	-----



### SHOULDER

- > Larger shoulder specific pad
- > Extended hose length
- > Ideal for use with torn Rotator Cuff, pinched nerve, Bursitis, Tendonitis, and contusions

			Single	20 Pack
Regular Hose	U	Non-Sterile	SV003006	SV103006



### HAND/WRIST/FINGERS

- > Universal left and right
- > Can be used for the wrist and hand or for the fingers

Regular Hose	U	Non-Sterile	SV003007
Regular Hose	U	Sterile	SV013007S






## COLD CONTROL PADS

### ANKLE

- > Universal left and right
- > Covers ankle and forefoot with one pad
- > Ideal for use with Achilles Tendonitis, fractures and sprains

			Single	20 Pack
Regular Hose	U	Non-Sterile	SV003003	SV103003



### 5 X 1 0 RECTANGLE

- > Can be used for a variety of applications
- > Angled hose for easy patient use
- > Ideal for use on the back

Regular Hose	U	Non-Sterile	SV003004
--------------	---	-------------	----------



### 1 0 X 1 4 RECTANGLE

- > Can be used for a variety of applications
- > Angled hose for easy patient use
- > Ideal for use on the back

Regular Hose	U	Non-Sterile	SV003005
--------------	---	-------------	----------



### SHOULDER WRAP

- > Can be used for a variety of applications
- > Angled hose for easy patient use
- > Use with Universal Pad or Shoulder Pad

U	Non-Sterile	SV006302
---	-------------	----------





## ALPHABETICAL INDEX

Part Number	Description	Suggested HCFA Code	Page	Part Number	Description	Suggested HCFA Code	Page
Custom	20.50 Patellofemoral Custom	L1846+L2800	17	RK0800	Crossover Front Closure, Neoprene		15
PF010	20.50 Patellofemoral OTS	L1845+L2795	17	RK0810	Crossover Front Closure, Neoprene w/Buttress		15
AL03601	Achilles Boot	L4386	37	RK0200	Crossover Pull-up, Neoprene		15
AL23601	Achilles Boot w/Air	L4360	37	RK0210	Crossover Pull-up, Neoprene w/Buttress		15
AL03511	Achilles Boot w/PTB	L4386	37	RK0900	Crossover ROM, Front Closure		23
AL03521	Achilles Boot w/PTB	L4386	37	RK0100	Crossover Short Pull-up, Neoprene		15
AL23511	Achilles Boot w/PTB and Air	L4360	37	RK0110	Crossover Short Pull-up, Neoprene w/Buttress		15
AL23521	Achilles Boot w/PTB and Air	L4360	37	RK0700	Crossover Short Front Closure, Neoprene		15
AL05300	AdjustaFit LC	L4386	31	RK0710	Crossover Short Front Closure, Neoprene w/Buttress		15
AL05600	AdjustaFit LC Mid-Calf	L4386	31	RK0630	Crossover Short, Front Closure		15
AL25600	AdjustaFit LC Mid-Calf w/Air	L4360	31	RK0640	Crossover Short, Pull-up		15
AL05700	AdjustaFit LC Mid-Calf w/Ankle-Heel Support	L4386	31	RK0650	Crossover Short, Pull-up w/Buttress		15
AL05400	AdjustaFit LC Pre-assembled	L4386	31	RK0663	Crossover Sport Hinges w/Straps		15
AL25400	AdjustaFit LC Pre-assembled w/Air	L4360	31	RK0600	Crossover, Front Closure		15
AL05405	AdjustaFit LC Pre-assembled w/Ankle-Heel	L4386	31	RK0610	Crossover, Pull-up		15
AL25300	AdjustaFit LC w/Air	L4360	31	RK0620	Crossover, Pull-up w/Buttress		15
AL05305	AdjustaFit LC w/Ankle-Heel Support	L4386	31	EK087000	Extender Plus	L1832	21
AL05360	AdjustaFit Wide	L4386	31	EK007	Extender Plus Shoe Insert		21
AL05660	AdjustaFit Wide Mid-Calf	L4386	31	EK088000	Extender Plus w/Cuffs	L1832	21
AL25660	AdjustaFit wide Mid-Calf w/Air	L4360	31	AL04300	FlatForm	L4386	37
AL05460	AdjustaFit Wide Pre-assembled	L4386	31	AL04310	FlatForm Mid-Calf	L4386	37
AL25460	AdjustaFit Wide Pre-assembled w/Air	L4360	31	AL24310	FlatForm Mid-Calf w/Air	L4360	37
AL25360	AdjustaFit Wide w/Air	L4360	31	AL24300	FlatForm w/Air	L4360	37
AR078	Aligner ESE Lateral	L1843	8	AL03602	Hi-Top w/Bootie	L4386	33
AR079	Aligner ESE Medial	L1843	8	AL23602	Hi-Top w/Bootie and Air	L4360	33
AE050100	ARC Shoulder		43	AL07202	Hi-Top w/Bootie and EZ Set Hinge	L4386	33
AE050300	ARC w/Pillow	L3670	44	AL27202	Hi-Top w/Bootie, EZ Set Hinge and Air	L4360	33
AE050200	ARC-XR Shoulder (External Rotation)		43	AL03600	Hi-Top w/Cuff	L4386	33
AS020	Axiom Ankle	L1971	27	AL23600	Hi-Top w/Cuff and Air	L4360	33
Custom	Axiom Custom	L1846	13	AL07200	Hi-Top w/Cuff and EZ Set Hinge	L4386	33
Custom	Axiom Custom, Magnesium	L1846	13	AL27200	Hi-Top w/Cuff, EZ Set Hinge and Air	L4360	33
Custom	Axiom OA	L1846	9	SA6090	Knee Support (Closed Patella)	L1825	15
PK114	Axiom OTS	L1845	13	SA6160	Knee Support (Open Patella)	L1825	15
PK119	Axiom OTS Athletic Version	L1845	13	SA611	Lateral Patella Stabilizer	L1800+L2795	18
MG119	Axiom OTS Athletic Version Magnesium	L1845	13	AK07000	Lever Lock	L1832	23
MG114	Axiom OTS Magnesium	L1845	13	AL08502	Lo-Top w/Bootie	L4386	33
Custom	Axiom-D Custom	L1846	12	AL28502	Lo-Top w/Bootie and Air	L4360	33
Custom	Axiom-D Custom Magnesium	L1846	12	AL07002	Lo-Top w/Bootie and EZ Set Hinge	L4386	33
PK118	Axiom-D OTS	L1845	12	AL27002	Lo-Top w/Bootie, EZ Set Hinge and Air	L4360	33
PK113	Axiom-D OTS Athletic Version	L1845	12	AL08500	Lo-Top w/Cuff	L4386	33
MG113	Axiom-D OTS Athletic Version Magnesium	L1845	12	AL28500	Lo-Top w/Cuff and Air	L4360	33
MG118	Axiom-D OTS Magnesium	L1845	12	AL07000	Lo-Top w/Cuff and EZ Set Hinge	L4386	33
AL05800	Bunion Boot	L4360	36	AL27000	Lo-Top w/Cuff, EZ Set Hinge and Air	L4360	33
AL25800	Bunion Boot w/Air Ankle Pad	L4360	36	ES040005	M-Brace Air	L1906	29
SA90100	Cast Shoe	L3260	40	ES042005	M-Brace Air (20 Pack)	L1906	29
AL042	Charcot Conformer	L4360	35	ES045005	M-Brace Air (50 Pack)	L1906	29
AL242	Charcot Conformer w/Air Ankle Pad	L4360	35	ES062005	M-Brace Combo (20 Pack, 10 Air, 10 Gel)	L1906	29
SA30100	Clavicle Support	L3650	45	ES065005	M-Brace Combo (50 Pack, 25 Air, 25 Gel)	L1906	29
SA20000	Clinic Arm Sling	A4565	44	ES050005	M-Brace Gel	L1906	29
SA20050	Clinic Shoulder Immobilizer	L3670	44	ES052005	M-Brace Gel (20 Pack)	L1906	29
AL041	Conformer	L4386	34	ES055005	M-Brace Gel (50 Pack)	L1906	29
AL241	Conformer w/Air	L4360	34	EK08000	Merit OR	L1832	22
AE055100	Cool ARC Shoulder		43	AL05500	Mid-Calf w/Cuff	L4386	33
AE055300	Cool ARC with Pillow	L3670	44	AL25500	Mid-Calf w/Cuff and Air	L4360	33
AE055200	Cool ARC-XR Shoulder (External Rotation)		43	AL07100	Mid-Calf w/Cuff and EZ Set Hinge	L4386	33
EK08300	Cool Merit OR	L1832	22	AL27100	Mid-Calf w/Cuff, EZ Set Hinge and Air	L4360	33



## ALPHABETICAL INDEX

Part Number	Description	Suggested HCFA Code	Page	Part Number	Description	Suggested HCFA Code	Page
AE012	Neutral Hand Position Attachment		47	TA018	Thruster RLF Lateral Short	L1843	7
SA9020	Night Splint	L4396	39	TA009	Thruster RLF Medial	L1843	7
AL04301	Nite N' Day	L4386	38	TA029	Thruster RLF Medial Athletic Version	L1843	7
AL24301	Nite N' Day w/Air	L4360	38	TA019	Thruster RLF Medial Short	L1843	7
Custom	OA Air, Custom	L1844	8	SA801	Trigger and Keeper Thumb Support	L3923	49
AJ005	Original Femoral Plate		22	Custom	Ultimate Ankle Custom	L1971	28
EK0250	Original Hinge and Thigh Attachment		22	AS13025	Ultimate Ankle, Left	L1971	28
AM005	Original Shoe Insert		22	AS23025	Ultimate Ankle, Right	L1971	28
AL005	Original Tibial Plate		22	Custom	Ultimate Dynamic Custom	L1846	14
EK0200	Original w/Standard Hinge	L1832	22	GK060	Ultimate Dynamic, ED	L1845	14
AL004	Original Winged Tibial Plate		22	JK220	Ultimate Dynamic, HA	L1845	14
SV003007	Pad Hand/Fingers w/Hand/Wrist Straps (Reg)		56	Custom	Ultimate Sport Custom	L1846	14
SV013007	Pad Hand/Fingers w/Hand/Wrist Straps (Reg), Sterile		56	Custom	Ultimate Sport Custom, HI	L1846	14
SV003003	Pad, Ankle w/Ankle Straps (Reg)		57	GK040	Ultimate Sport, ED	L1845	14
SV003005	Pad, Rectangle 10x14 w/Back Straps (Reg)		57	GK070	Ultimate Sport, ED, HI	L1845	14
SV003004	Pad, Rectangle 5x10 w/Univ Straps (Reg)		57	JK020	Ultimate Sport, HA	L1845	14
SV003106	Pad, Shoulder w/Univ Straps		56	JK120	Ultimate Sport, HA, HI	L1845	14
SV103106	Pad, Shoulder w/Univ Straps, 20 pack		56	SV005100	Unit, Cold Control		55
SV003102	Pad, Universal w/Univ Straps (EXT)		56	SV007106	Unit, Cold Control w/Shoulder Pad/Straps (EXT)		55
SV003002	Pad, Universal w/Univ Straps (Reg)		56	SV007102	Unit, Cold Control w/Univ Pad/Straps (EXT)		55
SV013002	Pad, Universal w/Univ Straps (Reg), Sterile		56	SV007002	Unit, Cold Control w/Univ Pad/Straps (Reg)		55
SV103002	Pad, Universal w/Univ Straps (Reg), 20 pack		56	SV107002	Unit, Cold Control w/Univ Pad/Straps (Reg), 6 pack		55
SV113002	Pad, Universal w/Univ Straps (Reg), 20 pack, Sterile		56	SV105100	Unit, Cold Control, 6 pack		55
SA9025	Padded PF Splint	L4397	39	SA201000	Universal Arm Slings	A4565	45
AT010	Pelvic Attachment		25	AB0140	Volare (Corset Only), 10 inch		53
ES007	PFS	L1902	38	AB0130	Volare (Corset Only), 8 inch		53
ES008	PFS (10 pack)	L1902	38	AB0080	Volare Chairback Frame		53
AT090	Philippin Post-Op Hip	L1686	25	AB0120	Volare Lite	L0626	53
AT093	Philippin Post-Op Package		25	AB00500	Volare Lordotic Pads		53
SA900	Post-op Shoe	L3260	40	AB0180	Volare LS, 10 inch	L0627	52
AL03200	ProGait ST	L4386	32	AB0170	Volare LS, 8 inch	L0627	52
AL05200	ProGait ST Mid-Calf	L4386	32	AB0230	Volare LS--XT, 10 inch	L0631	52
AL25200	ProGait ST Mid-Calf w/Air	L4360	32	AB0220	Volare LS--XT, 8 inch	L0631	52
AL05205	ProGait ST Mid-Calf w/Ankle-Heel Support	L4386	32	AB007004	Volare Panel, Mid, Ant/Post		53
AL07500	ProGait ST Mid-Calf w/EZ Set Hinge	L4386	32	AB007006	Volare Panel, XT, Posterior		53
AL27500	ProGait ST Mid-Calf w/EZ Set Hinge and Air	L4360	32	AB007002	Volare Panel, Low, Ant/Post		53
AL07505	ProGait ST Mid-Calf w/EZ Set Hinge and AHS	L4386	32	AB0280	Volare w/Chairback, 10 inch	L0637	51
AL23200	ProGait ST w/Air	L4360	32	AB0270	Volare w/Chairback, 8 inch	L0637	51
AL03205	ProGait ST w/Ankle-Heel Support	L4386	32	AL18500	Wee Walker	L4386	32
AL07400	ProGait ST w/EZ Set Hinge	L4386	32	SA7020	Wraprotor Ankle Stabilizer w/Speed Laces	L1902	29
AL27400	ProGait ST w/EZ Set Hinge and Air	L4360	32	SA7025	Wraprotor Ankle Stabilizer w/Standard Laces	L1902	29
AL07405	ProGait ST w/EZ Set Hinge and Ankle-Heel	L4386	32	SA802	Wrist Brace	L3908	48
SA6140	Reinforced Patella Stabilizer	L1800+L2795	18	SA803	Wrist Brace w/Thumb Splint	L3807	48
SV006302	Shoulder Wrap		57	SA800	Wrist Lacer	L3908	48
AS0024	Sidekick	L1971	28	SA805	Wrist Lacer With Thumb Spica	L3807	48
SA9030	Silicone Heel Cups (Pair)	L3485	40	Custom	Z-12 Brace, Custom	L1846	11
AT013	Simple Pelvic	L1686	25	AZ114	Z-12 Brace, OTS	L1845	11
SA6070	Single Panel Knee Immobilizer	L1830	24	AZ119	Z-12 EXT Brace, OTS	L1845	11
SA204000	Sling and Swathe Immobilizer	L3670	45	Custom	Z-12 OA Brace, Custom	L1846	11
AE012	Supination/Pronation Hand Attachment		47				
AE023	Telescoping Elbow	L3760	47				
SA60500	Tendon Strap	L3701(Elbow), L1825(Knee)	18 & 47				
SA6080	Three Panel Knee Immobilizer	L1830	24				
Custom	Thruster RLF Custom	L1844	7				
TA008	Thruster RLF Lateral	L1843	7				
TA028	Thruster RLF Lateral Athletic Version	L1843	7				



## PRODUCT NUMBER INDEX

Part Number	Description	Suggested HCFA Code	Page	Part Number	Description	Suggested HCFA Code	Page
Custom	20.50 Patellofemoral Custom	L1846+L2800	17	AL05305	AdjustaFit LC w/Ankle-Heel Support	L4386	31
Custom	Axiom Custom	L1846	13	AL05360	AdjustaFit Wide	L4386	31
Custom	Axiom Custom, Magnesium	L1846	13	AL05400	AdjustaFit LC Pre-assembled	L4386	31
Custom	Axiom OA	L1846	9	AL05405	AdjustaFit LC Pre-assembled w/Ankle-Heel	L4386	31
Custom	Axiom-D Custom	L1846	12	AL05460	AdjustaFit Wide Pre-assembled	L4386	31
Custom	Axiom-D Custom Magnesium	L1846	12	AL05500	Mid-Calf w/Cuff	L4386	33
Custom	OA Air, Custom	L1844	8	AL05600	AdjustaFit LC Mid-Calf	L4386	31
Custom	Thruster RLF Custom	L1844	7	AL05660	AdjustaFit Wide Mid-Calf	L4386	31
Custom	Ultimate Ankle Custom	L1971	28	AL05700	AdjustaFit LC Mid-Calf w/Ankle-Heel Support	L4386	31
Custom	Ultimate Dynamic Custom	L1846	14	AL05800	Bunion Boot	L4360	36
Custom	Ultimate Sport Custom	L1846	14	AL07000	Lo-Top w/Cuff and EZ Set Hinge	L4386	33
Custom	Ultimate Sport Custom, HI	L1846	14	AL07002	Lo-Top w/Bootie and EZ Set Hinge	L4386	33
Custom	Z-12 Brace, Custom	L1846	11	AL07100	Mid-Calf w/Cuff and EZ Set Hinge	L4386	33
Custom	Z-12 OA Brace, Custom	L1846	11	AL07200	Hi-Top w/Cuff and EZ Set Hinge	L4386	33
AB00500	Volare Lordotic Pads		53	AL07202	Hi-Top w/Bootie and EZ Set Hinge	L4386	33
AB007002	Volare Panel, Low, Ant/Post		53	AL07400	ProGait ST w/EZ Set Hinge	L4386	32
AB007004	Volare Panel, Mid, Ant/Post		53	AL07405	ProGait ST w/EZ Set Hinge and Ankle-Heel	L4386	32
AB007006	Volare Panel, XT, Posterior		53	AL07500	ProGait ST Mid-Calf w/EZ Set Hinge	L4386	32
AB0080	Volare Chairback Frame		53	AL07505	ProGait ST Mid-Calf w/EZ Set Hinge and AHS	L4386	32
AB0120	Volare Lite	L0626	53	AL08500	Lo-Top w/Cuff	L4386	33
AB0130	Volare (Corset Only), 8 inch		53	AL08502	Lo-Top w/Bootie	L4386	33
AB0140	Volare (Corset Only), 10 inch		53	AL18500	Wee Walker	L4386	32
AB0170	Volare LS, 8 inch	L0627	52	AL23200	ProGait ST w/Air	L4360	32
AB0180	Volare LS, 10 inch	L0627	52	AL23511	Achilles Boot w/PTB and Air	L4360	37
AB0220	Volare LS--XT, 8 inch	L0631	52	AL23521	Achilles Boot w/PTB and Air	L4360	37
AB0230	Volare LS--XT, 10 inch	L0631	52	AL23600	Hi-Top w/Cuff and Air	L4360	33
AB0270	Volare w/Chairback, 8 inch	L0637	51	AL23601	Achilles Boot w/Air	L4360	37
AB0280	Volare w/Chairback, 10 inch	L0637	51	AL23602	Hi-Top w/Bootie and Air	L4360	33
AE012	Neutral Hand Position Attachment		47	AL241	Conformer w/Air	L4360	34
AE012	Supination/Pronation Hand Attachment		47	AL242	Charcot Conformer w/Air Ankle Pad	L4360	35
AE023	Telescoping Elbow	L3760	47	AL24300	Flatform w/Air	L4360	37
AE050100	ARC Shoulder		43	AL24301	Nite N' Day w/Air	L4360	38
AE050200	ARC-XR Shoulder (External Rotation)		43	AL24310	Flatform Mid-Calf w/Air	L4360	37
AE050300	ARC w/Pillow	L3670	44	AL25200	ProGait ST Mid-Calf w/Air	L4360	32
AE055100	Cool ARC Shoulder		43	AL25300	AdjustaFit LC w/Air	L4360	31
AE055200	Cool ARC-XR Shoulder (External Rotation)		43	AL25360	AdjustaFit Wide w/Air	L4360	31
AE055300	Cool ARC with Pillow	L3670	44	AL25400	AdjustaFit LC Pre-assembled w/Air	L4360	31
AJ005	Original Femoral Plate		22	AL25460	AdjustaFit Wide Pre-assembled w/Air	L4360	31
AK07000	Lever Lock	L1832	23	AL25500	Mid-Calf w/Cuff and Air	L4360	33
AL004	Original Winged Tibial Plate		22	AL25600	AdjustaFit LC Mid-Calf w/Air	L4360	31
AL005	Original Tibial Plate		22	AL25660	AdjustaFit wide Mid-Calf w/Air	L4360	31
AL03200	ProGait ST	L4386	32	AL25800	Bunion Boot w/Air Ankle Pad	L4360	36
AL03205	ProGait ST w/Ankle-Heel Support	L4386	32	AL27000	Lo-Top w/Cuff, EZ Set Hinge and Air	L4360	33
AL03511	Achilles Boot w/PTB	L4386	37	AL27002	Lo-Top w/Bootie, EZ Set Hinge and Air	L4360	33
AL03521	Achilles Boot w/PTB	L4386	37	AL27100	Mid-Calf w/Cuff, EZ Set Hinge and Air	L4360	33
AL03600	Hi-Top w/Cuff	L4386	33	AL27200	Hi-Top w/Cuff, EZ Set Hinge and Air	L4360	33
AL03601	Achilles Boot	L4386	37	AL27202	Hi-Top w/Bootie, EZ Set Hinge and Air	L4360	33
AL03602	Hi-Top w/Bootie	L4386	33	AL27400	ProGait ST w/EZ Set Hinge and Air	L4360	32
AL041	Conformer	L4386	34	AL27500	ProGait ST Mid-Calf w/EZ Set Hinge and Air	L4360	32
AL042	Charcot Conformer	L4360	35	AL28500	Lo-Top w/Cuff and Air	L4360	33
AL04300	FlatForm	L4386	37	AL28502	Lo-Top w/Bootie and Air	L4360	33
AL04301	Nite N' Day	L4386	38	AM005	Original Shoe Insert		22
AL04310	Flatform Mid-Calf	L4386	37	AR078	Aligner ESE Lateral	L1843	8
AL05200	ProGait ST Mid-Calf	L4386	32	AR079	Aligner ESE Medial	L1843	8
AL05205	ProGait ST Mid-Calf w/Ankle-Heel Support	L4386	32	AS0024	Sidekick	L1971	28
AL05300	AdjustaFit LC	L4386	31	AS020	Axiom Ankle	L1971	27



# PRODUCT NUMBER INDEX

Part Number	Description	Suggested HCFA Code	Page	Part Number	Description	Suggested HCFA Code	Page
AS13025	Ultimate Ankle, Left	L1971	28	SA20000	Clinic Arm Slings	A4565	44
AS23025	Ultimate Ankle, Right	L1971	28	SA20050	Clinic Shoulder Immobilizer	L3670	44
AT010	Pelvic Attachment		25	SA201000	Universal Arm Slings	A4565	45
AT013	Simple Pelvic	L1686	25	SA204000	Slings and Swathe Immobilizer	L3670	45
AT090	Philippon Post-Op Hip	L1686	25	SA30100	Clavicle Support	L3650	45
AT093	Philippon Post-Op Package		25	SA60500	Tendon Strap	L3701 <sup>(Elbow)</sup> , L1825 <sup>(Knee)</sup>	18 & 47
AZ114	Z-12 Brace, OTS	L1845	11	SA6070	Single Panel Knee Immobilizer	L1830	24
AZ119	Z-12 EXT Brace, OTS	L1845	11	SA6080	Three Panel Knee Immobilizer	L1830	24
EK007	Extender Plus Shoe Insert		21	SA6090	Knee Support (Closed Patella)	L1825	15
EK0200	Original w/Standard Hinge	L1832	22	SA611	Lateral Patella Stabilizer	L1800+L2795	18
EK0250	Original Hinge and Thigh Attachment		22	SA6140	Reinforced Patella Stabilizer	L1800+L2795	18
EK08000	Merit OR	L1832	22	SA6160	Knee Support (Open Patella)	L1825	15
EK08300	Cool Merit OR	L1832	22	SA7020	Wraprator Ankle Stabilizer w/Speed Lacers	L1902	29
EK087000	Extender Plus	L1832	21	SA7025	Wraprator Ankle Stabilizer w/Standard Laces	L1902	29
EK088000	Extender Plus w/Cuffs	L1832	21	SA800	Wrist Lacer	L3908	48
ES007	PFS	L1902	38	SA801	Trigger and Keeper Thumb Support	L3923	49
ES008	PFS (10 pack)	L1902	38	SA802	Wrist Brace	L3908	48
ES040005	M-Brace Air	L1906	29	SA803	Wrist Brace w/Thumb Splint	L3807	48
ES042005	M-Brace Air (20 Pack)	L1906	29	SA805	Wrist Lacer With Thumb Spica	L3807	48
ES045005	M-Brace Air (50 Pack)	L1906	29	SA900	Post-op Shoe	L3260	40
ES050005	M-Brace Gel	L1906	29	SA90100	Cast Shoe	L3260	40
ES052005	M-Brace Gel (20 Pack)	L1906	29	SA9020	Night Splint	L4396	39
ES055005	M-Brace Gel (50 Pack)	L1906	29	SA9025	Padded PF Splint	L4397	39
ES062005	M-Brace Combo (20 Pack, 10 Air, 10 Gel)	L1906	29	SA9030	Silicone Heel Cups (Pair)	L3485	40
ES065005	M-Brace Combo (50 Pack, 25 Air, 25 Gel)	L1906	29	SV003002	Pad, Universal w/Univ Straps (Reg)		56
GK040	Ultimate Sport, ED	L1845	14	SV003003	Pad, Ankle w/Ankle Straps (Reg)		57
GK060	Ultimate Dynamic	L1845	14	SV003004	Pad, Rectangle 5x10 w/Univ Straps (Reg)		57
GK070	Ultimate Sport, ED, HI	L1845	14	SV003005	Pad, Rectangle 10x14 w/Back Straps (Reg)		57
JK020	Ultimate Sport, HA	L1845	14	SV003007	Pad Hand/Fingers w/Hand/Wrist Straps (Reg)		56
JK120	Ultimate Sport, HA, HI	L1845	14	SV003102	Pad, Universal w/Univ Straps (EXT)		56
JK220	Ultimate Dynamic, HA	L1845	14	SV003106	Pad, Shoulder w/Univ Straps (EXT)		56
MG113	Axiom-D OTS Athletic Version Magnesium	L1845	12	SV005100	Unit, Cold Control		55
MG114	Axiom OTS Magnesium	L1845	13	SV006302	Shoulder Wrap		57
MG118	Axiom-D OTS Magnesium	L1845	12	SV007002	Unit, Cold Control w/Univ Pad/Straps (Reg)		55
MG119	Axiom OTS Athletic Version Magnesium	L1845	13	SV007102	Unit, Cold Control w/Univ Pad/Straps (EXT)		55
PF010	20.50 Patellofemoral OTS	L1845+L2795	17	SV007106	Unit, Cold Control w/Shoulder Pad/Straps (EXT)		55
PK113	Axiom-D OTS Athletic Version	L1845	12	SV013002	Pad, Universal w/Univ Straps (Reg), Sterile		56
PK114	Axiom OTS	L1845	13	SV013007	Pad Hand/Fingers w/Hand/Wrist Straps (Reg), Sterile		56
PK118	Axiom-D OTS	L1845	12	SV103002	Pad, Universal w/Univ Straps (Reg), 20 pack		56
PK119	Axiom OTS Athletic Version	L1845	13	SV103106	Pad, Shoulder w/Univ Straps (EXT), 20 pack		56
RK0100	Crossover Short Pull-up, Neoprene		15	SV105100	Unit, Cold Control, 6 pack		55
RK0110	Crossover Short Pull-up, Neoprene w/Buttress		15	SV107002	Unit, Cold Control w/Univ Pad/Straps (Reg), 6 pack		55
RK0200	Crossover Pull-up, Neoprene		15	SV113002	Pad, Universal w/Univ Straps (Reg), 20 pack, Sterile		56
RK0210	Crossover Pull-up, Neoprene w/Buttress		15	TA008	Thruster RLF Lateral	L1843	7
RK0600	Crossover, Front Closure		15	TA009	Thruster RLF Medial	L1843	7
RK0610	Crossover, Pull-up		15	TA018	Thruster RLF Lateral Short	L1843	7
RK0620	Crossover, Pull-up w/Buttress		15	TA019	Thruster RLF Medial Short	L1843	7
RK0630	Crossover Short, Front Closure		15	TA028	Thruster RLF Lateral Athletic Version	L1843	7
RK0640	Crossover Short, Pull-up		15	TA029	Thruster RLF Medial Athletic Version	L1843	7
RK0650	Crossover Short, Pull-up w/Buttress		15				
RK0663	Crossover Sport Hinges w/Straps		15				
RK0700	Crossover Short Front Closure, Neoprene		15				
RK0710	Crossover Short Front Closure, Neoprene w/Buttress		15				
RK0800	Crossover Front Closure, Neoprene		15				
RK0810	Crossover Front Closure, Neoprene w/Buttress		15				
RK0900	Crossover ROM, Front Closure		23				





## WARRANTY

---

### Bledsoe Warranty

All Bledsoe products are warranted to be free from defects in material and workmanship for a period of 90 days for wraps, pads and normal wear components. Additional warranties include the following:

#### Custom Rigid Bracing:

Products warranted to be free from defects in material and workmanship for:

- > 6 months for pads, straps and normal wear components
- > 12 months on all other parts
- > For the life of the device against hinge or shell breakage

#### OTS Rigid Bracing:

Products are warranted to be free from defects in material and workmanship for:

- > 6 months for pads, straps and normal wear components
- > 12 months on all other parts
- > 5 years against hinge or shell breakage (*3 years for Z-12 braces*)

These warranties apply to devices that have not been modified or subjected to misuse, abuse or neglect.

## ORDERING INFORMATION

---

Customer Care: 1.888.BLED SOE (253.3763)

Hours: 7:00 am - 7:00 pm Central Time

Fax: 972.660.5495

orders@bledsoebrace.com

www.BledsoeBrace.com





2009™  
**Bledsoe**

WWW.BLEDISOEBRACE.COM





BRACING FOR THE FUTURE

.....

2601 Pinewood Drive • Grand Prairie • Texas • 75051  
Local Phone 972.647.0884 • Local Fax 972.660.5495  
1.888.BLED SOE (253.3763) • [www.BledsoeBrace.com](http://www.BledsoeBrace.com)



**Bledsoe**